

# Miele washer-disinfectors for CSSDs and hospitals





# Quality from Miele Professional – Never be satisfied with less!



In developing reprocessing solutions for hospitals and CSSDs, Miele Professional has attached great significance over the past 50 years to working closely both with specialists in the medical field and instrument manufacturers. A cooperation based on trust and spawning innovations which prove their value in workaday life.

# **50 YEARS Medical Technology**

#### Quality

Your patients rely on your professional skills, so why not in turn place your trust in Miele's proverbial quality! The combination of choice materials, perfect workmanship and mature and sophisticated technology has earned Miele its exceptional reputation in the field of medical technology. And Miele's after-sales service is renowned for its award-winning excellence.

#### Capacity

Machines from Miele Professional excel in terms of their superior performance. A wide selection of accessories and processes reflecting everyday, practical needs leads to perfect results and high throughput, even in highly specialised applications. This ensures that expensive instruments are ready for use again in next to no time.

#### **Efficiency**

Profit on a daily basis from the unique economic proposition machines from Miele Professional represent: They are thorough and fast and are also gentle on materials. Thanks to their high quality, they are also low-maintenance and durable, adding to the value of your investments. The sparing use of water and chemicals further contributes towards keeping operating costs in check.



Developed for daily use: System4Med, a new approach to professional instrument reprocessing, meets the highest of demands by addressing all aspects of the instrument reprocessing cycle. This systematic approach ensures perfect results whilst reducing the strain on materials and ensuring economic reprocessing in modern-day CSSDs.



- Washing and disinfection
   Innovative washer-disinfectors with customisable features.
- Sterilisation
   Powerful large steam sterilisers for the fast reprocessing of instruments and textiles.
- Documentation
   Intelligent documentation software for reproducible and uninterrupted process monitoring
- Guarantees:
   Miele guarantees the highest of standards in providing advice,
   financing, after-sales service and validation of the highest Miele
   quality.



# Customised systems in three size categories:







Miele stands for innovation with respect to systematic washing, disinfection and the systematic sterilisation of medical devices: All system components, processes, after-sales service – for the flexible, reliable and efficient use of validatable machine-based processes.

Please do not hesitate to contact us if you require further information:

- Planning negotiations
- Advisory services / Expertise
- Conceptualisation
- Economic feasibility calculations
- Installation plans

### Exclusive to Miele

This logo highlights products with special properties only available from Miele.

#### **Explanation of symbols**

AE Stainless-steel casing
AP Drain pump
AV Dump valve
BO Boiler
D Steam heating
D/EL Steam/Electric Convertible
EL Electric heating

GS Glass door
KD After-sales service
ORTHO OrthoVario programme
OXI OxiVario programme
PT Profitronic controls
TA Drying unit

#### **Contents**

| Cleaning programmes                     | Page 6  |
|---|---------|
| Service                                 | Page 7  |
| Washer-disinfectors PG 8527 and PG 8528 | Page 8  |
| Load carriers and inserts               | Page 24 |
| Transport trolleys                      | Page 44 |
| Washer-disinfectors G 7823 and G 7824   | Page 46 |
| Load carriers and inserts               | Page 54 |
| Washer-disinfectors G 7825 and G 7826   | Page 60 |
| Load carriers and inserts               | Page 68 |

Cf. separate brochure:

Miele's new container and trolley washers PG 8821, PG 8822, PG 8825, PG 8830

for the reprocessing of transport trolleys, sterile supply containers, theatre shoes and beds and the DGD 8803/8805 for mattress disinfection.

## Cleaning programmes

Miele innovations for particularly critical cases. Milestones in the optimisation of cleaning results:

#### 1994



The **Vario TD** method can now be seen as the standard programme for routine instrument cleaning and disinfection, achieving good removal of protein-based contamination (blood, secretion) from instruments which do not pose a particular challenge. Thermal disinfection takes place at > 90°C for 5 minutes. To protect instruments, the final rinse is performed ideally with demineralised water without surfactant. The Vario TD-AN programme with its modified parameters is available to reprocess thermo-stable anaesthetic instruments.

- Intensive cleaning below protein denaturation temperature
- Disinfection according to EN ISO 15883
- Exceptional material compatibility

#### 2004



#### OxiVario

Special programme for the critical instruments used in trauma surgery as well as high-frequency cauterising instruments requiring higher standards of cleaning.

- Excellent cleaning and removal of organic soil
- Time-saving as there is no need for preand post-treatment



#### OxiVario Plus

Special programme to prevent the iatrogenic spread of vCJD according to guidelines published by the task force set up by Germany's Robert Koch Institute.

- Excellent cleaning and removal of organic soil
- Time-saving as there is no need for preand post-treatment

#### 2005



#### OrthoVario

Special programme for orthopaedic instruments including motor systems and other medical products made from aluminium.

- Excellent cleaning performance
- Good material compatibility even on instruments sensitive to alkalinity

#### 2011



#### RobotVario

Special-purpose programmes, matching process chemicals and load carriers for the reprocessing of instruments frequently used in robotic-assisted surgery.

#### Since 2014: RobotVario X-tra

From the end of 2014 onwards, the newly developed ROBOTVARIO X-tra system will further shorten cycle times and also offers increased capacity and allows the reprocessing of various generations of robotic-assisted instruments.

- Accommodates up to 8 shafted instruments and offers space for two further special instruments
- Modular design of load carrier enables reprocessing of various instrument generations and types
- Thorough irrigation of all lumens thanks to innovative injector system
- Integrated microfilter for additional filtration of wash liquor

### Service



#### Highly skilled Miele after-sales service

Miele offers service packages which are closely geared to the requirements of CSSDs.

A team of skilled and knowledgeable medical product advisers is available to support clients with their expertise both during the planning phase and later during operation.

In addition to a wide range of services, Miele service engineers also provide assistance in optimising and redesigning cleaning processes and reprocessing methods. With its product lineup, Miele always tries to contribute towards simplifying procedures in hospitals and making them safer and more efficient.

#### Services, Service contracts

In order to guarantee maximum forward cost visibility, Miele recommends that its clients take out service contracts:

- Inspection
- Maintenance
- Full service maintenance

#### Validated performance

Miele also offers a series of process checks performed by highly qualified Miele medical product engineers in accordance with the relevant legislation and standards as well as national guidelines and good-practice recommendations:

- Initial validation consists of installation, operation and performance qualification after machine installation
- Revalidation (repeat performance qualification) is generally necessary every 12 months, after maintenance, repairs or after any modifications to installation or operation parameters.

#### Exclusive to Miele

- Blanket service coverage combined with fast response
- 90% of service calls result in first-time fixes
- Key genuine spares available for 15 years after discontinuation of series production

It is not without reason that Miele's after-sales service operation has been acclaimed for many years in succession for its excellence (in an annual survey performed by ServiceBarometer AG, Munich).





# Washer-disinfectors PG 8527 and PG 8528


# Adding value to instrument reprocessing in CSSD

Miele Professional has been setting milestones for decades in the field of efficient and safe machine-based reprocessing in clinics and surgeries. And, once again, Miele is setting new standards with the PG 85 machine generation. When it comes to centrally or decentrally reprocessing large quantities of instruments, Miele's PG 8527 and PG 8528 offer huge benefits: Greater cleaning capacity, improved process security and better efficiency.



| Performance comparison<br>Washer-disinfectors | Unit width/depth | Doors                 | Cabinet<br>Useable dimensions<br>H/W/D<br>Useable volume | Approx. capacity per cycle                   |
|---|------------------|-----------------------|--|--|
| PG 8527 / PG 8528                             | 1150/870 mm      | Vertical-rise door(s) | 675/650/800 mm<br>351 l                                  | 18 DIN mesh trays<br>7 AN sets<br>3 MI sets  |
| G 7823 / G 7824                               | 900/770 mm       | Bottom-hinged door(s) | 510/530/620 mm<br>168 l                                  | 8 DIN mesh trays<br>3 AN sets:<br>2 MI sets  |
| G 7825 / G 7826                               | 900/750 mm       | Bottom-hinged door(s) | 683/541/610 mm<br>225 l                                  | 10 DIN mesh trays<br>4 AN sets:<br>2 MI sets |

# PerfectEco - combined water and energy savings



# PG85 Perfect Eco

- Economic operation through lower water and electricity consumption
- Excellent washing and disinfection

Thanks to an innovative combination of water and heat recovery, large-chamber decontamination units from the PG 85 series can be put to use in CSSDs even more efficiently. Waste heat from the drying cycle is channeled through the steam condenser to pre-heat the incoming demineralised water supply. The preheated demineralised water is then heated to a higher temperature in the boiler in readiness for the final rinse. This way, process heat is kept within the system and put to direct use. A heat-recovery steam condenser also dispenses with the need for an intake of cold water which, on the unconverted model, is discharged after use into the drainage system. So in addition to energy savings, this approach results in a significant reduction in water consumption of up to 130 I per cycle.

PerfectEco is also available as a retrofit kit.

- 1 Demineralised water intake, cold
- 2 Heat-recovery steam condenser
- 3 Demineralised water boiler with additional heating
- 4 Demineralised water intake into chamber
- 5 Hot waste air
- 6 Cooled waste air

# PerfectTouchControl – Simple, local-language controls PerfectPure sensor – Residue-free rinsing

## PG85 | Perfect Touch Control

- Optimum user convenience
- Reliable hygiene
- Perfect control

Simple to operate, easy to clean: Washer-disinfectors from the PG 85 series feature a touch-sensitive display. This easy-to-use PerfectTouch display guarantees unique user convenience combined with superb hygiene. A fully flush, chemical-proof display screen makes for simple and effective wipe disinfection.

The controls are outlined on the glass surface and slight pressure is enough to activate functions and launch programmes, even when wearing protective gloves.

The man-machine interface involves the use of only a very limited number of controls; all steps in the process appear in the display in the user's own language. Display texts, for example for A0 values, actual temperatures, conductivity, countdown times and all other protocol data can be defined individually. Using the data, integrated A0 value controls, parameters can be set to comply with national standards or to suit a hospital's own regime.



- Chemical-resistant glass surfaces
- Innovative programme cycles



### PG85 | PerfectPureSensor

- Continuous conductivity monitoring
- Residue-free rinse quality
- Safe reprocessing

Residue in the final rinse can have a negative impact on reprocessing results. Alkaline residue in ophthalmology, for instance, can in extreme cases result in serious complications when instruments are reused. Excessive residue can also result in corrosion and deposits on instruments.


The new and patented PerfectPure conductivity meter reliably detects the presence of minerals in the rinse water, such as the dissolved salts introduced with alkaline or acidic process chemicals, limiting them to a threshold level defined by the user.

Conductivity monitoring is contact-free using a maintenance-free system with extremely precise tolerance in the 5–40  $\mu$ S/cm and 40  $\mu$ S/cm - 100 mS/cm ranges. Depending on machine settings, conductivity readings can even be used to

control the programme cycle. This allows the number of necessary rinse cycles to be determined automatically if conductivity is outside the predetermined range: Recording conductivity over the entire process helps ensure the reproducibility of validated processes. Recording conductivity over the entire process helps ensure the reproducibility of validated processes.


### Exclusive to Miele

• Maintenance-free conductivity monitoring





# CombiEco System – The intelligent approach to water and energy savings in CSSD



With the help of a manifold, water from the final rinse on PG 8528 washer-disinfectors is fed into a buffer tank serving a container and trolley washer.

#### Saving resources with Miele

Miele washer-disinfectors from the PG 8528 series and container and trolley washers from the PG 88 series already achieve excellent standards of cleaning and disinfection combined with the efficient use of water and energy. Thanks to the new cross-product approach to water and energy recycling afforded by the CombiEco system, both can now be deployed even more efficiently in the CSSD context.

In a typical CSSD with several PG 8528 units and a container and trolley washer, valuable fully demineralised water can be collected after the rinse cycles on PG 8528 washer-disinfectors and fed via a manifold to a buffer tank to enable the heat energy in the water to be used in the container and trolley washer when used to reprocess for example containers and transport trolleys. This saves on valuable and expensive resources and further considerably improves a CSSD's green credentials.

# **PerfectFlow sensor** – Volumetric monitoring of dispensed quantities **PerfectSpeed sensor** – Spray arm monitoring

### PG85 | PerfectFlowSensor

- Continuous monitoring of dispensed volumes
- Precise results, user-defined tolerances
- Perfect dispensing control

A decisive factor contributing to good reprocessing results is the precise volumetric control of dispensed chemicals. The currently valid DIN EN ISO 15883 standard prescribes the redundant monitoring of dispensing. Miele's new PerfectFlow sensor using ultrasound technology offers considerably greater safety margins than conventional systems. The PerfectFlow sensor is a standard feature and guarantees a hitherto unparallelled degree of precision in controlling and monitoring volumetric flow, independent of viscosity and ambient temperatures. The monitoring system is fully independent of the dispensing system and can be calibrated.

Dispensing tolerances can be set to comply with legal standards or defined individually; chemicals are dispensed efficiently and reliably, irrespective of the type of product or ambient conditions (continuous operation, fluctuating climatic conditions). Any deviation from the target quantities are safely detected and the reproducibility of validated processes guaranteed. An error message is issued or the programme is aborted if values are outside the tolerance range.



### Exclusive to Miele

- Extremely precise
- Product- and temperature-independent dispensing control

# PG85 | Perfect SpeedSensor

- Precisely monitored reprocessing
- Improved reproducibility of validated processes
- Perfect washing and disinfection results

To guarantee perfect and safe cleaning and disinfection results, the rotational speed of the spray arms must be within defined limits. With the new PerfectSpeed sensor, the precise speed of each individual spray arm is carefully monitored and documented - whether in the cabinet or on board baskets and load carriers. The spray arm monitoring feature uses a sensor strip located outside the cabinet to detect the passage of spray arms and to ensure that speeds are within the specified range. Information shown in the display indicates whether the values are correct or whether the user must intervene on account, for example, of excessive foam slowing spray arm motion down.

In the event of a deviation from target values, either an error message is issued or the programme is interrupted immediately to allow the user to deal with the cause of the fault, depending on system parameters.

Deviations can also be recorded by automatic process documentation. Spray arm sensing offers effective protection against spray arm blockages by items in the load and also provides information on pressure conditions in the machine, on load carriers and in baskets. And, most importantly, spray arm speeds provide a valuable indication as to the reproducibility of validated processes, increasing safety margins in machine-based instrument reprocessing systems by a considerable degree.



### Exclusive to Miele

- Spray arm monitoring on all levels
- Monitoring of spray arms speeds

# PerfectFineFilter – Micro-filtration of wash liquor PerfectDoc – Gap-free documentation

## PG85 | Perfect FineFilter

- Micro-fine filtration with 0.2 mm mesh size
- · Highly efficient filtration of wash liquor

The thorough removal of residuals from the wash liquor equates to greater security in reprocessing, better instrument functionality and the preservation of the value of instruments.

With the micro-fine filter with a mesh size of only 0.2 mm combined with an expansive surface area, even the smallest of particles are reliably removed from circulation.

Particularly when treating the narrow-lumened instruments used in ophthalmology and micro-surgery, users stand to profit greatly from better cleaning performance and a reduced need for inspection.



## Exclusive to Miele

 Micro-fine filtration with 0.2 mm mesh size

# PG85 | Perfect Doc

- Continuous process documentation
- Wide range of parameters documented
- Perfect process traceability

The PG 8527/28 is fitted with a network interface for process documentation as a standard feature. The PerfectDoc module allows the machine to interface with process documentation software. This facilitates the recording of many process parameters, including temperatures, as well as the documentation of entire process protocols including A0 values, dispensed quantities, spray arms speeds and conductivity readings. Alternatively, a printer can be used for documentation purposes.



### Exclusive to Miele

- Simple integration into standard process documentation systems without additional hardware or software
- Comprehensive sterile supplies logistics with process documentation, instrument management and administration of sterile supplies

## Machine features and specifications



#### Versions

- PG 8527: Single-door model with vertical sliding door
- PG 8528: Two-door barrier washerdisinfector for installation in a diaphragm wall

#### Capacity per cycle

 18 DIN mesh trays or 6 DIN containers 150 x 300 x 600 mm or 4 DIN containers 300 x 300 x 600 mm or 7 AN sets or 3 MI sets

#### Design

- Stand-alone or side-by-side installation
- Width 1150 mm
- Modular machine concept, customisable features
- Single-chamber system for washing, disinfection and drying
- Service-friendly design
- Heater elements outside wash chamber
- Low heat and sound emissions thanks to double insulation

# Vertical-rise all-glass doors and cabinet lighting

To improve process visibility and supervision, all-glass vertical-rise doors are available as optional extras on the new range of products. The door can be programmed to open automatically at the end of a programme. Door motion is monitored electronically.

Models with all-glass doors also feature cabinet lighting comprising 4 integrated spotlights. Lighting times are programmable.

#### Cleaning technology

- Hygienic freshwater system with fresh water intake in each programme stage
- Cleaning, disinfection and drying in a closed, single-cabinet system
- Hygienic stainless-steel wash cabinet with coved corners and sloping selfdraining ceiling
- 2 spray arms in cabinet for thorough cleaning of instrument surfaces
- · Spray arms with high water impact force
- Full water jet access, ensuring optimum results

- Thorough cleaning of lumens thanks to injector system
- Direct docking of load carriers to water circuit

#### Standard features

- 2 powerful circulation pumps
- Triple filtration with large surface filter, coarse filter and micro-fine filter
- Filters in inlet hoses
- Flowmeter to monitor water intake quantities
- 1 dump valve
- Spray arm monitoring
- Ultrasound flow control
- Automatic load carrier detection

#### Dispenser systems

- 2 bellows-type dispenser pumps for liquid detergent and neutralising agent
- Extendible (further 3 internal dispenser pumps retrofittable)



• Maintenance-free dispensing technology

Technical data on Pages 22/23



#### Water softeners

• High-capacity water softener (external)

#### Drainage

• 2 drain pumps

#### **Controls**

- Freely programmable controls PROFITRONIC
- 64 programme slots
- 16 standard washing and disinfection programmes
- 17 service programmes
- 31 vacant programme slots
- User interface with local-language display
- Display of programme selection and programming dialogs, countdown time, faults and operating hours.
- Compilation of new programmes using machine controls or using PC/laptop via optical interface

#### Interfaces

- 1 x RJ 45 network interface for process documentation
- 4 serial interfaces, for example for remote service
- Optical interface for service and maintenance

#### **Printer**

• Integrated 8-pin printer to record key process parameters

#### **Automatic load carrier recognition**

Sensors for load carrier recognition and programme assignment

#### Multiport

• For connection of printer and/or scanner

#### Safety features

- Electric door lock
- Programme recontinuation in event of power outage
- Peak-load negotiation
- Optical and acoustic signal at end of programme
- 2 separate sensors for temperature monitoring and control
- Port for simple positioning of sensors in the wash cabinet for process validation
- Sensors in cabinet and a magnetic strip on load carriers for automatic assignment of programmes to load carriers
- Volume flow control
- Spray arm sensing

# Washer-disinfectors



# Machine versions, additional modules


| PG 8527  | Features                | Mat. no. | Art. no.   |
|--|-------------------------|----------|------------|
| Electric   | Features<br>AE PT EL AV | 6881680  | 62.8527.21 |
| Steam/Electric   | AE PT D/EL AV           | 6881690  | 62.8527.21 |
| Steam/Electric   | AE PT D/EL AV           | 0001090  | 02.0027.31 |
| PG 8528  |                         |          |            |
| Electric   | AE PT EL AV             | 6757350  | 62.8528.21 |
| Steam/Electric   | AE PT D/EL AV           | 6757370  | 62.8528.31 |
| Explanation of symbols cf. Page 5  |                         |          |            |
| Additional modules for PG 8527 and PG 8528   |                         | Mat. no. | Art. no.   |
| Drain pump   |                         | 6758120  | 69.2400.01 |
| Cabinet AISI 316 L (DIN 1.4404) grade stainless steel  |                         | 6758160  | 69.2410.01 |
| Boiler EL  |                         | 6758220  | 69.2430.01 |
| Boiler D/EL  |                         | 6758270  | 69.2430.02 |
| Glass door for PG 8527 (1 door)  |                         | 6758320  | 69.2450.01 |
| Glass doors for PG 8528 (2 glass doors)  |                         | 6815810  | 69.2450.02 |
| Printer – integrated, on infeed or outfeed side (please quote)   |                         | 6758340  | 69.2470.01 |
| Scanner connection with scanner on infeed, outfeed or both sides (please quote)  | 1                       | 7686510  | 69.2470.12 |
| Conductivity module (information on Page 41)   |                         | 6758400  | 69.2440.01 |
| Dispenser pump – Integrated (additional) <sup>2</sup> – Please state purpose (cf. footnotes)   |                         | 6758410  | 69.2460.01 |
| OxiVario   |                         | 6758420  | 69.2460.20 |
|  |                         |          |            |
| Please note:  D/EL machine version can only be operated in combination with D/EL boiler  Combination of printer/scanner only possible on outfeed side  | hemical disinfectant.   | 6815800  | 69.2460.30 |
| Please note:  D/EL machine version can only be operated in combination with D/EL boiler  Combination of printer/scanner only possible on outfeed side  | hemical disinfectant.   | 6815800  | 69.2460.30 |
| Please note:  D/EL machine version can only be operated in combination with D/EL boiler  Combination of printer/scanner only possible on outfeed side  | nemical disinfectant.   | 6815800  | 69.2460.30 |
| Please note:  D/EL machine version can only be operated in combination with D/EL boiler  Combination of printer/scanner only possible on outfeed side  | hemical disinfectant.   | 6815800  | 69.2460.30 |
| Please note:  D/EL machine version can only be operated in combination with D/EL boiler  Combination of printer/scanner only possible on outfeed side  | nemical disinfectant.   | 6815800  | 69.2460.30 |
| Please note:  D/EL machine version can only be operated in combination with D/EL boiler  Combination of printer/scanner only possible on outfeed side  | hemical disinfectant.   | 6815800  | 69.2460.30 |
| Please note:  D/EL machine version can only be operated in combination with D/EL boiler  Combination of printer/scanner only possible on outfeed side  | nemical disinfectant.   | 6815800  | 69.2460.30 |
| Please note:  D/EL machine version can only be operated in combination with D/EL boiler  Combination of printer/scanner only possible on outfeed side  | nemical disinfectant.   | 6815800  | 69.2460.30 |
| Please note:  D/EL machine version can only be operated in combination with D/EL boiler  Combination of printer/scanner only possible on outfeed side  | hemical disinfectant.   | 6815800  | 69.2460.30 |
| OrthoVario (programme package also contains OxiVario)  Please note:  D/EL machine version can only be operated in combination with D/EL boiler  1 Combination of printer/scanner only possible on outfeed side  2 Please specify whether dispenser pump is for detergent, neutralising agent or ch | nemical disinfectant.   | 6815800  | 69.2460.30 |
| Please note:  D/EL machine version can only be operated in combination with D/EL boiler  Combination of printer/scanner only possible on outfeed side  | hemical disinfectant.   | 6815800  | 69.2460.30 |
| Please note:  D/EL machine version can only be operated in combination with D/EL boiler  Combination of printer/scanner only possible on outfeed side  | nemical disinfectant.   | 6815800  | 69.2460.30 |
| Please note:  D/EL machine version can only be operated in combination with D/EL boiler  Combination of printer/scanner only possible on outfeed side  | nemical disinfectant.   | 6815800  | 69.2460.30 |
| Please note:  D/EL machine version can only be operated in combination with D/EL boiler  Combination of printer/scanner only possible on outfeed side  | hemical disinfectant.   | 6815800  | 69.2460.30 |
| Please note:  D/EL machine version can only be operated in combination with D/EL boiler  Combination of printer/scanner only possible on outfeed side  | nemical disinfectant.   | 6815800  | 69.2460.30 |
| Please note:  D/EL machine version can only be operated in combination with D/EL boiler  Combination of printer/scanner only possible on outfeed side  | hemical disinfectant.   | 6815800  | 69.2460.30 |
| Please note:  D/EL machine version can only be operated in combination with D/EL boiler  Combination of printer/scanner only possible on outfeed side  | nemical disinfectant.   | 6815800  | 69.2460.30 |
| Please note:  D/EL machine version can only be operated in combination with D/EL boiler  Combination of printer/scanner only possible on outfeed side  | nemical disinfectant.   | 6815800  | 69.2460.30 |

# Components/Accessories



### SBW

- Plinth/floor tray
   Frame with integrated stainlesssteel floor tray
  • 2 cross-brackets for moving
- machine
- Cut-outs for steam and plumbing connections, dump valve, dual drain valve and floor tray drain
- H 100, W 1150, D 856 mm



#### SBW/1 Plinth/drip tray without cutouts

- For version with drain pump
- Utilities, e.g. water, electricity and steam, through ceiling
   H 100, W 1150, D 856 mm



#### SBWR Plinth/floor tray

- Roller plinth/floor tray
- Castors allow machine to be pulled forward for servicing
- PG 8527 with drain pump
- Utility supply access through ceiling • H 100, W 1150, D 856 mm



#### Integrated printer for process documentation

- 8-needle printer and RS 232 serial interface
  • Paper width: 58 mm
- For installation on infeed side of PG 8527 or outfeed/infeed side of PG 8528
- of PG 8528
   Recording of following parameters during programme: Date and Mach. no., Prog. no., prog. name, starting and finishing time, dispensing concentration, dispensing temperature and pump 1–4, target temperature reached (wash/dry) with times, all faults (e.g. 'water inlet'), all incidents of manual intervention (start, stop, power outage) (start, stop, power outage)
- Meets process parameters with respect to temperature and holding time

Paper rolls available from Miele Spares: Mat. no.: 4781470

Paper roll size: Width 58 mm, Ø 50 mm, Length 20 m

| Accessories | Features   | Mat. no. | Art. no.   |
|-------------|--|----------|------------|
| SBW         | Plinth/drip tray for PG 8527/PG 8528                   | 6757850  | 69.2530.01 |
| SBW/1       | Plinth/drip tray for PG 8527/PG 8528, without cut-outs | 6757860  | 69.2530.02 |
| SBWR        | Plinth/drip tray, on castors for PG 8527               | 5653140  | 69.3710.05 |
| Printer     | For process documentation                              | 6758340  | 69.2470.01 |

### Accessories



### Drying unit/electric

- Side-channel compressor
- Drying of interior and exterior load surfaces
- 2 pre-filters Class EU 4, filter rating > 95% (ASHRAE Standard 52–68)
- Life cycle 200 h
- 2 particulate filters/HEPA H 13, Filter rating > 99.95% (DIN EN 1822-2011)
- Life cycle 1000 hVoltage 3N AC 400 V 50 Hz

- Heating 2 x 4 kW = 8 kW
  Heavy-duty fan/side-channel compressor, life cycle above 10,000 h, 2 air circuits 1.8 kW
- Total rated load 10 kW
- Air throughput approx. 250 m³/h
- Temperature infinitely adjustable between 60 and 115°C
- Time programmable between 1 and 240 mins.
- Incl. fitting frame for installation on top of PG 8527/PG 8528
- Panelling through to ceiling to be provided on site, cf. MAV 27/28



## Drying unit/steam

- Side-channel compressor
- Drying of interior and exterior load surfaces
- 2 pre-filters Class EU 4, filter rating > 95% (ASHRAE Standard 52-68)
- Life cycle 200 h
- 2 particulate filters/HEPA H 13, Filter rating > 99.95% (DIN EN 1822-2011)
- Life cycle 1000 h
- Heating steam-to-air heat exchanger (steam circuit stainless steel, air circuit aluminium)
- Steam pressure 3.5-6 bar (350-600 kPa)
- Steam quality: filtered saturated steam
- Steam capacity (max.) 15 kg/h (machine and drying unit 50 kg/h)
- Heavy-duty fan/side-channel compressor, life cycle above 10,000 h, 2 air circuits 1.8 kW
- Total rated load 1.8 kW
- Air throughput approx. 250 m<sup>3</sup>/h
- Temperature infinitely adjustable between 60 and 115°C
- Time programmable between 1 and 240 mins.
- Incl. fitting frame for installation on top of PG 8527/PG 8528
- Panelling through to ceiling to be provided on site, cf. MAV 27/28



#### DK 27/28 Heat-exchanger steam condenser

- Water-cooled (only dehumidified air should be introduced into air conditioning system)
- · Connection to on-site cooling circuit (no water consumption) or cold water connection (water consumption)
- Max. water pressure: 8 bar
- On-site installation
- · Reduction of air discharge temperature to approx. 30-35°C
- Reduction in relative humidity to approx. 60-70%

#### DK 27/28 WR

PerfectEco heat-recovery steam condenser

- Preheats demineralised water in steam condenser
- Water savings up to 130 I/cycle
- On-site installation
- Further information, cf. Page 10



#### **MAV 27/28** Fitting kit/top-box panelling kit for drying unit/steam condenser frame

- · Service hatches, with lock, for both sides of diaphragm wall, stainless steel
- Ventilation grille on the infeed side
- H 760, W 1150, D 765 mm
- Diaphragm wall panelling between the top of the top-box panelling and the ceiling must be
- provided on site

  Incl. MAV top panel to cover topbox panelling

| Accessories | Features  | Mat. no. | Art. no.   |
|-------------|---|----------|------------|
| TA/E        | Drying unit/electric                                  | 6757710  | 69.2500.01 |
| TA/D        | Drying unit/steam                                     | 6757770  | 69.2500.02 |
| DK 27/28    | Heat-exchanger steam condenser                        | 6757790  | 69.2510.01 |
| DK 27/28 WR | Heat-recovery steam condenser**                       | 9032080  | 69.2511.01 |
| MAV 27/28   | Top-box panelling kit                                 | 6757820  | 69.2520.01 |
| OxiVario    | OxiVario retrofit kit incl. pumps and installation    | 6911330  |            |
| OrthoVario  | OrthoVario retrofit kit incl. pumps and installation* | 6911320  |            |

<sup>\*</sup> presumes installation of OxiVario

<sup>\*\*</sup> only for systems with boiler

# Technical data

| Washer-disinfectors  | PG 8527 (single-door mode | l) PG 8528 (barrier model |
|--|---------------------------|---------------------------|
| Model with vertical sliding doors  | •                         | -                         |
| Through-feed unit with vertical-lift doors                                 | -                         | •                         |
| All-glass doors/cabinet lighting   | 0                         | 0                         |
| Single/multiple installations  | •                         | •                         |
| Freshwater system, max. temperature 93°C                                   | •                         | •                         |
| Direct load carrier docking for cleaning and drying of lumened instruments | •                         | •                         |
| 2 circulation pumps [Qmax. I/min]  | 400/600*                  | 400/600*                  |
| Boiler for heating demineralised water                                     | 0                         | 0                         |
| Controls/Programmes  |                           |                           |
| Profitronic +, 16 standard programmes                                      | •                         | •                         |
| 64 programme slots   | •                         | •                         |
| Electric door lock   | •                         | •                         |
| Peak-load negotiation  | •                         | •                         |
| Network interface for process documentation                                | •                         | •                         |
| Magnetic strip for automatic load carrier recognition                      | •                         | •                         |
| Spray arm sensing  | •                         | •                         |
| Conductivity metering  | 0                         | 0                         |
| Remote service enabled   | •                         | •                         |
| Water connections  |                           |                           |
| 1 x cold water, 2–10 bar flow pressure (200–1000 kPa) (max. 4°dH)          | •                         | •                         |
| 1 x hot water, 2–10 bar flow pressure (200–1000 kPa) (max. 4°dH)           | •                         | •                         |
| 1 x demineralised water, 2–10 bar flow pressure (200–1000 kPa)             | •                         | •                         |
| 3 inlet hoses ½" with ¾" threaded union                                    | •                         | •                         |
| Drain valve DN 50, odour trap to be fitted on site                         | •                         | •                         |
| 2 drain pumps, DN 22, odour trap to be provided on site                    | 0                         | 0                         |
| Electrical connection: Electric heating                                    |                           |                           |
| 3 N AC 400 V 50 Hz   | •                         | •                         |
| Cabinet heating [kW]   | 18                        | 18                        |
| Boiler heating [kW]  | 15                        | 15                        |
| Circulation pump [kW]  | 0.7/1.2*                  | 0.7/1.2*                  |
| Total rated load w/o TA [kW]   | 20                        | 20                        |
| Total rated load with electric TA [kW]                                     | 20                        | 20                        |
| Fuse rating [A]  | 3 x 32                    | 3 x 32                    |
| Electrical connection: Steam heating                                       |                           |                           |
| 3 N AC 400 V 50 Hz   | •                         | •                         |
| Circulation pump [kW]  | 0.7/1.2*                  | 0.7/1.2*                  |
| Total rated load w/o TA [kW]   | 2                         | 2                         |
| Total rated load with steam-heated TA [kW]                                 | 2                         | 2                         |
| Total rated load with electric TA [kW]                                     | 10                        | 10                        |
| Fuse rating [A]  | 3 x 16                    | 3 x 16                    |
| Steam connection G ½" (DN 15)  | •                         | •                         |
| Operating pressure 350-600 kPa (steam-heated drying unit)                  | •                         | •                         |
| Compressed air connection 600–1200 kPa                                     | •                         | •                         |
| Electrical connection: Steam/Electric (convertible)                        |                           |                           |
| 3 N AC 400 V 50 Hz   | •                         | •                         |
| Cabinet heating [kW]   | 18                        | 18                        |
| Boiler heating [kW]  | 15                        | 15                        |
| Circulation pump [kW]  | 0.7/1.2*                  | 0.7/1.2*                  |
| Total rated load with electric TA [kW]                                     | 20                        | 20                        |
| Fuse rating [A]  | 3 x 32                    | 3 x 32                    |
| Steam connection G ½" (DN 15)  | •                         | •                         |
| Operating pressure 350–1000 kPa (electric drying unit)                     | •                         | •                         |
| Compressed air connection 600–1200 kPa                                     | •                         | •                         |
|  |                           |                           |
|  |                           |                           |
| * Spray arms / direct docking  |                           |                           |

# Technical data

| Nasher-disinfectors  | PG 8527 (single-door model) | PG 8528 (barrier mo |
|--|-----------------------------|---------------------|
| Dispenser systems  |                             |                     |
| 2 bellows-type pumps for detergent and neutralising agent                                      | •                           | •                   |
| 2 x 10 l supply canisters  | •                           | •                   |
| Space for 4 x 10 l supply canisters  | •                           | •                   |
| /olume flow control  | •                           | •                   |
| Max. 3 additional dispenser pumps*   | 0                           | 0                   |
| Dimensions/Weight  |                             |                     |
| External dimensions H incl. plinth tray [mm]   | 1660                        | 1660                |
| External dimensions H incl. frame with TA [mm]   | 2420                        | 2420                |
| External dimensions W/D [mm]   | 1150/870                    | 1150/870            |
| Jseable cabinet dimensions H/W/D [mm]  | 675/650/800                 | 675/650/800         |
| Overall cabinet dimensions H/W/D [mm]  | 860/685/800                 | 860/685/800         |
| Loading height above floor   | 850                         | 850                 |
| Veight [kg]  | 408                         | 386                 |
| Outer casing   |                             |                     |
| Stainless steel (AE)   | •                           | •                   |
| Test certificates  |                             |                     |
| /DE, VDE-EMC, IP 20, MDD CE 0366   | •                           | •                   |
| A/E drying unit, electrically heated   |                             |                     |
| Supply voltage   | 3 N AC 400 V 50 Hz          | 3 N AC 400 V 50 Hz  |
| an [kW]  | 1.8                         | 1.8                 |
| Heater bank, depending on model [kW]   | 8                           | 8                   |
| otal rated load, depending on model [kW]   | 10                          | 10                  |
| Air throughput [m³/h]  | approx. 250                 | approx. 250         |
| [emperature selection in 1° increments [°C]  | 60–115                      | 60–115              |
| Fime selection in 1-minute increments  | 1–240                       | 1–240               |
| 2 X pre-filters EU 4, filter rating > 95%, filter life 200 h                                   | •                           | •                   |
| 2 particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-2011), filter life 1000 h | •                           | •                   |
| TA/D drying units steam-heated   |                             |                     |
| Steam pressure [bar]   | 3.5–6                       | 3.5-6               |
| Steam quality: filtered saturated steam  | •                           | •                   |
| Steam output: (max.) [kg/h]  | 15                          | 15                  |
| machine and drying unit 50 kg/h)   |                             |                     |
| Total rated load, depending on model [kW]  | 1.8                         | 1.8                 |
| Air throughput [m³/h]  | approx. 250                 | approx. 250         |
| [emperature selection in 1° increments [°C]  | 60–115                      | 60–115              |
| Fime selection in 1-minute increments  | 1–240                       | 1–240               |
| 2 x pre-filters EU 4, filter rating > 95%, filter life 200 h                                   | •                           | •                   |
| 2 particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-2011), filter life 1000 h | •                           | •                   |
|  |                             |                     |
| 4 pumps required for OxiVario programme, cf. additional module on Page 21.                     |                             |                     |
| 5 pumps required for OxiVario programme, cf. additional module on Page 21.                     |                             |                     |
|  |                             |                     |
|  |                             |                     |
|  |                             |                     |
|  |                             |                     |
|  |                             |                     |
|  |                             |                     |
|  |                             |                     |
|  |                             |                     |
|  |                             |                     |
|  |                             |                     |
|  |                             |                     |
|  |                             |                     |
|  |                             |                     |
|  |                             |                     |
|  |                             |                     |
|  |                             |                     |

### Load carriers with 2 - 6 levels



#### E 975/2 load carrier (empty)

- For inserts on 2 levels
- Load capacity:6 DIN mesh trays
- Built-in spray arm
- Loading dimensions (from bottom):
   Level 1: H 297, W 590, D 780 mm
   Level 2: H 290, W 590, D 780 mm
- Magnetic strip for automatic load carrier recognition (incl. magnets)
- H 427, W 640, D 790 mm

Mat. no. 7 765 790 Art. no. 69.5975.03



#### E 900-5/2 load carrier (empty)

- For inserts on 5 levels
- Load capacity:15 DIN mesh trays
- 4 built-in spray arms
- Clearances from bottom:
   Level 1: H 80, W 585, D 780 mm
   Levels 2–4: H 80, W 595, D 780 mm
   Level 5: H 73, W 595, D 780 mm
- Magnetic strip for automatic load carrier recognition (incl. magnets)
- H 605, W 640, D 790 mm

Mat. no. 7 765 760 Art. no. 69.5900.06



#### E 935/2 load carrier (empty)

- For inserts on 3 levels
- Load capacity:9 DIN mesh trays
- 2 built-in spray arms
- Clearances from bottom:
   Level 1: H 202, W 585, D 780 mm
   Level 2: H 202, W 595, D 780 mm
   Level 3: H 132, W 595, D 780 mm
- Magnetic strip for automatic load carrier recognition (incl. magnets)
- H 524, W 640, D 790 mm

Mat. no. 7 765 780 Art. no. 69.5935.03



#### E 900-6 load carrier (empty)

- For inserts on 6 levels, 2 thereof removable
- Connections for lumened instruments
- Load capacity:18 DIN mesh trays
- 3 built-in spray arms
- Clearances from bottom:
   Levels 1, 3: H 71 (162), W 555, D 780 mm
   Levels 2 and 4: H 71, W 565, D 780 mm
   Level 5: H 71, W 580, D 780 mm
   Level 6: H 71, W 595, D 780 mm
- Magnetic strip for automatic load carrier recognition (incl. magnets)
- H 553, W 636, D 790 mm

Mat. no. 9 525 120 Art. no. 69.6900.07



#### E 900-4/2 load carrier (empty)

- For inserts on 4 levels
- Load capacity:12 DIN mesh trays
- 3 built-in spray arms
- Clearances from bottom:
   Level 1: H 112.5, W 585, D 780 mm
   Level 2: H 112.5, W 595, D 780 mm
   Level 3: H 112.5, W 595, D 780 mm
   Level 4: H 114, W 595, D 780 mm
- Magnetic strip for automatic load carrier recognition (incl. magnets)
- H 557, W 640, D 790 mm

Mat. no. 7 765 740 Art. no. 69.5900.05

### Exclusive

- Spray arm monitoring on all levels
- Monitoring of spray arms speeds
- Max. load per mesh tray 10 kg
- Offers highest capacity, removable levels and connections for lumened instruments on one and the same load carrier

# Load carrier with 4 levels and injectors



# E 901/2 load carrier with drying connection (empty)

- For mesh trays, bowls and 8 lumened instruments
- 3 built-in spray arms
- 4 levels
- Clearances from bottom:
   Level 1: H 80, W 593, D 780 mm
   Level 2: H 80, W 593, D 780 mm
   Level 3: H 80, W 593, D 780 mm
   Level 4: H 210, W 520, D 780 mm
- Load capacity for levels 1-3:
   3 DIN mesh trays (485 x 254 x 50 mm)
   or
- 3 DIN mesh trays (540 x 254 x 50 mm) or
- 2 ISO mesh trays (485 x 344 x 50 mm)
- Load capacity Level 4:
- 3 DIN mesh trays (485 x 254 x 50 mm) or
- 2 DIN mesh trays (540 x 254 x 50 mm) or
- various bowls, injector nozzles for lumened instruments
- Connection for hot-air drying unit
- Magnetic strip for automatic load carrier recognition (incl. magnets)
- H 590, W 640, D 790 mm

#### Supplied as standard with:

- 3 x E 362 blanking screws
- 2 x E 447 adapters, female, for male LuerLock
- 2 x E 449 adapters, male, for LuerLock, female
- 4 x injection nozzles 4 mm, length 110 mm (Mat. no. 4 175 030)
- 3 x injection nozzles 4 mm, length 120 mm with spring support (Mat. no. 4 398 951)
- 3 x E 980 angled adapters

Mat. no. 6 893 870 Art. no. 69.5901.03



#### E 980 angled adapter

- Enables connection of long lumened instruments
- Screw fittings for various LuerLock adapters or nozzles for lumened instruments

Mat. no. 7 025 360 Art. no. 69.5980.01


### Inserts for theatre shoes



#### E 984 insert 1/2

- For various utensils
- For holders
- -7 x E 985 for 12 kidney dishes or
- 5 x E 986 for 5 bowls or
- 5 x E 987 for 15 theatre shoes max. Size 44 or
- 7 x E 989 universal holders for 8 insoles or
- -3 x E 993 for 9 theatre shoes, max. Size 48
- H 65 (150), W 371, D 572 mm

Mat. no. 5 337 750 Art. no. 69.5984.01



#### E 984 insert with E 987

• E 984 can be fitted with 5 x E 987 holders, each for 15 theatre shoes, max. Size 44

#### E 987

Mat. no. 5 337 800 Art. no. 69.5987.01



#### E 984 insert with E 989

• E 984 can be fitted with 7 x E 989 universal holders, e.g. for 8 insoles

#### E 989

Mat. no. 5 376 740 Art. no. 69.5989.01



#### E 930 insert for theatre shoes

- For E 975/2 load carrier on lower level
- 36 holders, 280 mm for 36 theatre shoes, max. Size 43
- H 289, W 586, D 752 mm

Mat. no. 5 001 360 Art. no. 69.5930.01



#### E 931 insert for theatre shoes

- For E 975/2 load carrier on lower and upper levels
- 20 holders for 20 theatre shoes max. Size 48
- Total capacity per cycle:
   40 theatre shoes (20 pairs)
- H 254, W 586, D 752 mm

Mat. no. 7 439 490 Art. no. 69.5931.01

## Inserts for kidney dishes, bowls Test Kit



E 492 insert 1/2

- For 9 kidney dishes
- 9 holders, H 86 mm, spacing 49 mm
- H 120, W 256, D 474 mm


Mat. no. 5 771 140 Art. no. 69.5492.01



E 984 insert with E 985

Can be fitted with 7 holders
 E 985 for 12 kidney dishes

Mat. no. 5 337 760 Art. no. 69.5985.01



E 984 insert with E 986

• Can be fitted with 5 holders E 986 for 5 bowls

Mat. no. 5 337 770 Art. no. 69.5986.01



E 936 insert for boom lamp handles

- For various lamp handles
- 15 receptacles
- L 500, W 245, H 90 mm

Mat. no. 7 759 480 Art. no. 69.5936.01



#### TK/1 test kit

- Tests for the presence of proteins and monitors cleaning results
- Contents sufficient for 48 tests
- With coding strips for reflectometer

#### Post-reprocessing security

In collaboration with Merck, Miele has developed a protein test (test kit) to provide a fast method of checking cleaning performance on instruments. This allows regular quality control of cleaning performance.

Mat. no. 6 157 330



#### ProCare protein check

Test to determine protein residue after thermal disinfection

- Proof of protein residue and qualitative evaluation
- Precise application with indicator pen
- Short 3-minute exposure time
- 3 pens per box

### Load carrier for containers



E 910/4 load carrier (empty)

- For 6 DIN containers and lids
- Container size: 150 x 300 x 600 mm
- Magnetic strip for automatic load carrier recognition (incl. magnets)
- H 455, W 640, D 790 mm

Mat. no. 7 439 500 Art. no. 69.5910.05



#### E 911/3 load carrier (empty)

- For 4 DIN containers and lids
- Container size: 300 x 300 x 600 mm
- Magnetic strip for automatic load carrier recognition (incl. magnets)
- H 455, W 640, D 790 mm

Mat. no. 6 943 850 Art. no. 69.5911.04



E 912/2 load carrier (empty)

- For 4 ISO containers and lids
- Container size: 180 x 400 x 600 mm
- Magnetic strip for automatic load carrier recognition (incl. magnets)
- H 455, W 640, D 790 mm

Mat. no. 6 943 860 Art. no. 69.5912.03



#### E 913/3 load carrier with drying connection (empty)

- Load carrier for containers/lids/mesh trays 3 levels
- Load dimensions: Level 1: for 3 DIN containers 150 x 300 x 600 mm with lid
- Level 2: H 97, W 575, D 780 mm Level 3: H 103, W 595, D 780 mm
- 2 built-in spray arms
- Connections for 5 lumened instruments
- Connection for hot-air drying unit
- Magnetic strip for automatic load carrier recognition (incl. magnets)
- H 574, W 640, D 790 mm

Mat. no. 6 893 880 Art. no. 69.5913.04



E 934/1 load carrier (empty)

- For orthopaedic containers on 4 levels
- Clearances from bottom: Level 1: H 90 or 180 mm (w/o Level 2), W 585, D 780 mm Level 2: H 90, W 595, D 780 mm Level 3: H 167, W 595, D 780 mm Level 4: H 180, W 595, D 780 mm
- Level 2 is removable. This gives a clearance on Level 1 of 180 mm.
- Spray arms below Levels 3 and 4
- Magnetic strip for automatic load carrier recognition (incl. magnets)
- H 500, W 640, D 790 mm

Mat. no. 7 765 770 Art. no. 69.5934.02

# Load carriers and modules for orthopaedic instruments



# E 941 load carrier with drying connector (empty)

- For modules on 2 levels
- Depending on the size of the glassware, up to 2 modules can be accommodated on each level
- Water and drying air enters via a direct docking system and adapters.
- Loading dimensions (from bottom upwards):

Level 1 = (without top module) H 609, W 558, D 352 mm Level 1 = (with top module) H 317, W 558, D 352 mm Level 2

H 245, W 558, D 352 mm

- Connection for hot-air drying unit
- Magnetic strip for automatic load carrier recognition (incl. magnets)
- H 421, W 619, D 790 mm

Mat. no. 4 812 530 Art. no. 69.5941.01



Miele ORTHOVARIO approved for the reprocessing of current-series Aesculap motor systems.



#### E 948 injector module for E 941

- For orthopaedic instruments
- For E 142 mesh tray (not included)
- Equipped with 28 nozzles
- Loading dimensions:
   H 125, W 520, D 310 mm
- H 230, W 558, D 352 mm

Mat. no. 7 025 310 Art. no. 69.5948.01

#### E 945/2 module (not illustrated)

- Carriage frame for inserts
- H 55, W 558, D 352 mm

Mat. no. 7 459 420 Art. no. 69.5945.03



#### E 949 injector module for E 941

- For marrow-channel pins and reamers
- 8 sleeves for marrow channel pins (prefitted) and 6 extension sleeves for bone marrow drills (supplied)
- H 115, W 558, D 352 mm

Mat. no. 7 025 320 Art. no. 69.5949.01



#### E 955 filter insert

- For lumened instruments, e.g. sleeves with diameter of 3-15 mm and max.
   length of 125 mm
- H 85, W 221, D 156 mm

Mat. no. 7 025 330 Art. no. 69.5955.01

### Load carriers for micro-instruments



Geuder AG recommends that their current range of instruments be reprocessed in a Miele washer-disinfector.





# E 929/1 mobile injector unit with drying connector

- 3 levels for approx. 6 OP sets
- Vertical clearance: Level 1: 94 mm
   Level 2: 100 mm
   Level 3: 287 mm
- Capacity: 6 DIN mesh trays
- 45 nozzles on the right-hand side (23 male and 22 female LuerLock adapters)
- 28 nozzles on left-hand side (horizontal), each with hose and male LuerLock adapter (E 792)
- Connection for 2 x E 478 on feed pipe
- Connection for hot-air drying unit
- Magnetic strip for automatic load carrier recognition (incl. magnets)
- H 492, W 640, D 790 mm

Mat. no. 6 893 890 Art. no. 69.5929.02

#### Supplied as standard with:

- 14 connectors LuerLock female/female (E 790)
- 14 connectors LuerLock female/male (E 791)
- 2 x E 476 adapters
- 1 x E 477 stopper

#### E 478/1 adapter

• For 4 narrow-lumened cannulae (Sautter cannulae)

Mat. no. 7 292 920 Art. no. 69.7478.02



#### FP ED

- Stainless-steel filter insert for E 478/1
- Ø 30 mm

Mat. no. 7 290 770 Art. no. 69.7512.02

# Load carriers for anaesthetic equipment



# E 915/3 mobile injector unit with drying connector

- For approx. 7 AN sets or 20 breathing tubes, max. length 1.5 m, and various accessories
- Tubes wound around spiral drum
- Direct injection system
- Connection for hot-air drying unit
- Magnetic strip for automatic load carrier recognition
- H 637, W 640, D 790 mm

#### Supplied as standard with:

- 1 x E 142 insert 1/2 mesh tray, H 45 (55), W 255, D 480 mm
- 10 x E 466 injector nozzles for breathing bags, 8 x 333 mm
- 18 x E 496 injector nozzles, 4 x 120 mm
- 1 x E 916 holder for 5 corrugated breathing tubes (supplied)
- 4 x E 917 holders for 5 silicone breathing tubes
- 1 x E 918 holder for 5 paediatric breathing tubes (supplied)
- 1 x E 431 injector nozzle for breathing bags, (supplied)
- 20 spring-loaded 12 mm supports for breathing tubes, max. 1.5 m
- 1 magnetic strip for automatic load carrier recognition

Mat. no. 6 943 870 Art. no. 69.5915.04



#### E 497 adapter

- Nozzle, small
- For double- and single-tubed laryngeal masks
- H 195 mm, with side pipe, 109 mm
- Illustration on right Mat. no. 7 759 470 Art. no. 69.7497.01

#### E 498 adapter

- Nozzle, large
- For double- and single-tubed laryngeal masks
- H 238 mm, with side pipe, 120 mm
- Illustration on left

Mat. no. 7 760 120 Art. no. 69.7498.01



#### E 466 accessory for AN load carrier

Injector nozzle for breathing bag, 8 x 333 mm

Mat. no. 4 966 620 Art. no. 69.7466.01

# E 496 accessory for AN load carrier (no illustration)

Injector nozzle for intubation material,
 4 x 120 mm

Mat. no. 5 235 490 Art. no. 69.7496.01

# Load carrier for MI instruments and reprocessing dental transmission instruments




E 909 load carrier with drying connector (empty)

- Modular system for approx. 1-2 MI sets
- Levels 1–3 vacant for mesh trays,
  e.g. E 142
- Clearances from bottom:
   Level 1: H 80, W 585, D 780 mm
   Level 2: H 80, W 595, D 780 mm
   Level 3: H 80, W 595, D 780 mm
   Level 4: H 216, W 515, D 780 mm
- Connection for hot-air drying unit
- Magnetic strip for automatic load carrier recognition (incl. magnets)
- H 698, W 640, D 790 mm

Mat. no. 7 201 500 Art. no. 69.5909.01

### Exclusive to Miele

 Simultaneous cleaning of mesh trays and MI OP sets



#### E 925 module

- For short MI instruments
- For up to 20 lumened instruments, max. length 190 mm in top row and max. length 220 mm in bottom row
- H 40, W 568, D 510 mm

Mat. no. 7904100 Art. no. 69.5925.01



#### E 926 module

- For long shafted MI instruments
- For up to 13 lumened instruments
- H 40, W 568, D 510 mm

Mat. no. 7 297 050 Art. no. 69.5926.01

#### Supplied as standard with E 926:

- 1 x E 336 MIBO irrigation sleeve, Length 121 mm, 11 mm
- 2 x E 362 blanking screws, 11 mm
- 5 x E 442 irrigation sleeves, length 121 mm, 11 mm, with 6 mm cap
- 3 x E 443 irrigation sleeves, length 121 mm, 11 mm, with 10 mm cap
- 2 x E 448 silicone tubes with LuerLock adapter, female
- 1 x E 454 injector nozzles for trocar sleeves, 10-15 mm
- 1 x E 464 sleeves for E 454 injector nozzles
- 2 x E 456 spring stays for MI instruments

#### Supplied as standard with E 925:

- 1 x E 336 MIBO irrigation sleeve, Length 121 mm, 11 mm
- 2 x E 362 blanking screws, 11 mm
- 1 x E 442 irrigation sleeve, length 121 mm, 11 mm, with cap 6 mm
- 2 x E 447 female adapter
- 4 x E 448 hose with thread
- 2 x E 449 LuerLock adapter, male
- 4 x E 452 injector nozzles, 2.5 x 60 mm, screw thread
- 3 x E 453 injector nozzles, 4.0 x 110 mm with clamps, screw thread
- 3 x E 454 injector nozzles for trocar sleeves, 10-15 mm, screw thread
- 1 x E 464 adapter for injector nozzle E 454
- 1 x E 472 clamp for injector nozzle, Ø 4 mm
- 1 x E 907 mesh basket with lid, for sundry small items, with lid

# Load carrier for MI instruments and reprocessing dental transmission instruments



# E 902/1 mobile injector unit with drying connector

- Modular system for approx. 3 MI OP sets
- Accommodates up to 78 lumened instruments in 6 modules with integrated nozzles/adapters.
- Top and bottom levels free for mesh trays, e. g. E 142 Upper level: H 64, W 595, D 780 mm Lower level:
- H 53, W 585, D 780 mm
   Connection for hot-air drying unit
- Magnetic strip for automatic load carrier recognition (incl. magnets)
- H 620, W 640, D 790 mm

#### Supplied as standard with:

- 3 x E 905/1 modular insert for short minimally invasive instruments
- 3 x E 906/1 modular inserts for long minimally invasive instruments
- 3 x E 444 insert/spiral rack for fibre optic cables or tubing
- 2 x E 445, 12 caps for E 442 irrigation sleeves

6 mm orifice

• 1 x E 446, 12 caps for E 443 irrigation sleeves,

10 mm orifice

- 2 m silicone hose, Ø 5 mm Mat. no. 4 822 830
- 6 x E 362 blanking screws

Mat. no. 6 059 530 Art. no. 69.5902.02

### Exclusive to Miele

 Simultaneous cleaning of mesh trays and MI OP sets

# Modules for E 902/1 load carrier



#### E 903/1 modular insert

- For TUR sets (transurethral resection)
- For up to 10 lumened instruments
- H 40, W 461, D 510 mm

#### Supplied as standard with:

- 3 x E 442 irrigation sleeves, 121 mm, for MI instruments, 4–8 mm
- 1 x E 444 insert/spiral rack for fibre optic cables and irrigation tubes
- 1 x E 447 adapter, female for LuerLock, male
- 4 x E 448 silicone hoses, 300 mm, 5 x 1.5 mm
- 3 x E 453 injector nozzles,4.0 x 110 mm with spring supports
- 1 x E 454 injector nozzles for trocar sleeves, 10-15 mm
- 3 x E 467 injector sleeves, 205 mm, for MI instruments/forceps
- 3 x E 469 injector sleeves, 300 mm, for MI instruments/urology
- 1 x E 907/1 insert / mesh tray with lid for small parts
- 2 m silicone hose, 5 mm Mat. no. 4 822 830
- 1 plastic support

Mat. no. 6 943 230 Art. no. 69.5903.02



#### E 905/1 modular insert

- For short MI instruments
- For up to 16 lumened instruments
- Sectioned for arthroscopes, laparoscopes etc.
- H 40, W 461, D 510 mm

#### Supplied as standard with:

- 1 x E 336 MIBO irrigation sleeves for pipettes/MI instruments,
   11 mm, Length 121 mm
- 2 x E 362 blanking screws
- 1 x E 442 irrigation sleeve for MI instruments,

Length 121 mm, 11 mm, with 6 mm cap

- 2 x E 447 adapters, female adapter for LuerLock, male
- 4 x E 448 silicone tube with LuerLock adapter, female
- 2 x E 449 male adapter, for female LuerLock
- 4 x E 452 injector nozzles, 2.5 x 60 mm
- 3 x E 453 injector nozzles,
  4.0 x 110 mm with spring supports
- 3 x E 454 injector nozzles for trocar sleeves, 10-15 mm
- 1 x E 464 sleeve for E 454 injector nozzles
- 1 x E 907/1 insert / mesh tray with lid for small parts
- 1 x E 472 clamp for injector nozzle, 4 mm

Mat. no. 6 943 240 Art. no. 69.5905.02



#### E 906/1 modular insert

- For long shafted MI instruments
- For up to 10 lumened instruments
- Sectioned for arthroscopes, laparoscopes etc.
- H 40, W 461, D 510 mm

#### Supplied as standard with:

- 2 x E 362 blanking screws
- 5 x E 442 irrigation sleeves for MI instruments
- 3 x E 443 irrigation sleeves for MI instruments, 8–12 mm
- 2 x E 448 silicone tubes with LuerLock adapter
- 1 x E 454 injector nozzles for trocar sleeves, 10-15 mm
- 2 x E 456 spring stays for MI instruments
- 1 x E 464 sleeve for E 454 injector nozzles
- 1 x E 908/1 insert for separable MI instruments/inner shafts
- 1 x E 336 irrigation sleeve MIBO

Mat. no. 6 943 250 Art. no. 69.5906.02

#### Note

The modules can also be used on module units E 574 (for G 7823/G 7824) and E 774/1 (for G 7825/G 7826).

## Module for E 902/3 for reprocessing dental instruments



Leading manufacturers recommend the reprocessing of transmission instruments in Miele thermal disinfectors.



KaVo. Dental Excellence.



#### to Miele

#### E 919 dental module

- For connection of 10 turbines/handpieces and contra-angles as well as dental and ENT instruments
- Can be used on E 902/3 (new) as well as E 902/1, E 774/1, E 574, E 474/4
- Items should be reprocessed in the VARIO TD programme
- Cleaning should be performed using neutral to mildly alkaline liquid detergents
- Combination with MI instruments possible

Note: 10 AUF supplied. ADS must be ordered separately.

Mat. no. 7 759 490 Art. no. 69.5919.01

#### E 902/3 load carrier (empty)

Mat. no. 6 943 830 Art. no. 69.5902.04

#### Note

The module can also be used on E 574 (for G 7823/G 7824) and E 774/1 (for G 7825/G 7826).



#### AUF 1 for G 7881 and G 7891

- Sleeve for transmission instruments in O 177/1 upper basket
- Consists of: Seal, threaded bush, bracket, filter housing – lower section, filter housing – upper section
- 5 filter inserts

#### **AUF 2 for G 7831**

- Sleeve for transmission instruments in O 801/2 upper basket
- Consists of: Seal, threaded bush, bracket, filter housing – lower section, filter housing – upper section
- 5 filter inserts



FP filter inserts for AUF 1 and AUF 2

- 20 off
- Porosity 2
- •Ø 30 mm

# For connection of transmission instruments

#### for AUF 1 and AUF 2

- ADS 1 (white): Ø approx. 20 mm
- ADS 2 (green): Ø approx. 16 mm
- ADS 3 (red): Ø approx. 22 mm

## ROBOTVARIO X-tra -

## Centralised reprocessing of robotic-assisted instruments



A 502 load carrier with A 602 module and PG 8528 decontamination unit

#### Reprocessing with PG 8528

As a systematic approach to reprocessing robotic-assisted instruments, ROBOTVARIO enjoys great trust and confidence in specialist circles and is recommended by Intuitive Surgical Inc., manufacturer of robotic-assisted systems. Miele Professional has further optimised the process and now offers an increased load capacity and a modular load carrier system.

The new A 502 load carrier permits the use of the new ROBOTVARIO X-tra system together with various modules in centralised instrument reprocessing in hospitals. The load carrier is able to take up to 8 shafted instruments and 2 further special instruments, hold them securely and guarantee the thorough internal irrigation of all lumens. An additional row of nozzles ensures the targeted cleaning of the heavily soiled distal tips of instruments, particularly the fine crevices.

Furthermore, the A 502 load carrier is also equipped with a micro-filter which prevents particulate deposits collecting inside instruments. Alongside the shafted and special instruments, the same load carrier can also accommodate further instruments such as the other utensils used in robotic-assisted surgery, e.g. trocar sleeves.

#### A 502 load carrier

- For use on PG 8528
- Reprocessing of up to 8 shafted instruments used in robotic-assisted surgery and 2 special instruments
- Holders for further instruments and accessories
- Micro-filter with 8 cm diameter
- H 466, W 640, D 790 mm Art. no. 69.6502.01 Mat. no. 09530400

For use in A 502 load carrier:

#### A 600 module (not illustrated)

- For instruments used in robot-assisted surgery, older generations (S / Si)
- H 303, W 555, D 318 mm Mat. no. 10 131 910 Art. no. 69.6600.01

#### A 602 module

- For instruments used in robot-assisted surgery, most recent generation (Xi)
- H 303, W 555, D 318 mm Mat. no. 10 131 930 Art. no. 69.6602.01

#### A 805 connector (not illustrated)


- For special instruments, surgical staplers (Si)
- H 92, W 41, D 165 mm Mat. no. 10 131 950 Art. no. 69.6805.01

#### A 806 connector (not illustrated)

- For special instruments with side flush port
- H 77, W 41 mm, D 165 mm Mat. no. 10 131 960 Art. no. 69.6806.01

#### A 807 connector (not illustrated)

- For special instruments, forceps, older generation (S / Si)
- H 77, W 41, D 165 mm Mat. no. 10 131 980 Art. no. 69.6807.01



The A 502 load carrier features a modular design, a sophisticated injector system and a microfilter for the thorough, particle-free cleaning of the lumens on shafted instruments.

A 502, fitted with A 602 module



A 502 load carrier



Replacement of micro-filter in injector system



A 602 module



Accommodates robotic-assisted instruments

### Cross-product accessories



#### E 469 irrigation sleeve

- For MI instruments/urology
- Length 300 mm, 11 mm
- Cap, 6 mm
   (Mat. no. 4174960)
- Clamp (Mat. no. 4 174 850)

Mat. no. 5 563 950 Art. no. 69.5469.01

#### E 467 irrigation sleeve

- For MI instruments/forceps
- Length 205 mm, 11 mm Mat. no. 5 039 640 Art. no. 69.5467.01

#### E 336 injector sleeve MIBO

• Length 121 mm, 11 mm

Mat. no. 3 809 390 Art. no. 69.7336.01

#### E 442 irrigation sleeve

- For MI instruments, 4-8 mm, screw thread
- Length 121 mm, 11 mm
- Cap, 6 mm orifice (Mat. no. 4 174 960)
- Clamp (Mat. no. 4 174 850)

Mat. no. 4 203 830 Art. no. 69.7442.01

#### E 443 irrigation sleeve

- For MI instruments, 8-8.5 mm
- Length 121 mm, 11 mm
- Cap, 10 mm orifice (Mat. no. 4 174 970)
- Clamp (Mat. no. 4 174 850)

Mat. no. 4 203 840 Art. no. 69,7443.01



#### E 445 caps

- 12 caps for irrigation sleeves
- 6 mm orifice Mat. no. 4 203 850 Art. no. 69.7445.01

#### E 446 caps

- 12 caps for irrigation sleeves
- 10 mm orifice Mat. no. 4 203 860 Art. no. 69.7446.01



#### E 456 spring stay

• For MI instruments Mat. no. 4 203 960 Art. no. 69.7456.01

#### E 475 holder bracket

- For MI insert
- To stabilise irrigation sleeves (Mat. no. 4 692 430)

Mat. no. 5 226 040 Art. no. 69.7475.01



#### E 447 adapter, female

 For LuerLock, male, screw thread for E 450/1, O 176, O 183
 Mat. no. 4 203 870
 Art. no. 69.7447.01

### E 517 adapter, male, with centre pin

 For LuerLock, female, screw thread for E 450/1, O 176, O 183
 Mat. no. 7 790 040
 Art. no. 69.7517.01

## E 449 adapter, male, without centre pin

 For LuerLock, female, screw thread for E 450/1, O 176, O 183
 Mat. no. 4 203 890
 Art. no. 69.7449.01

#### E 452 injector nozzle

 Length 60 mm, Ø 2.5 mm, screw-thread, for mobile injector unit
 Mat. no. 4 203 920

Art. no. 69.7452.01

#### E 453 injector nozzle

- Clip bracket
- Length 110 mm, 4 mm, screw-thread, for mobile injector unit

Mat. no. 4 203 930 Art. no. 69.7453.01

#### E 448 silicone hose

- LuerLock adapter, male
- Length 300 mm, 5 mm
- Nozzle, 8 x 1 mm thread Mat. no. 4 203 880 Art. no. 69.7448.01



#### E 464 adapter

- For injector nozzle E 454,
  13 x 65 mm
- Height-adjustment clip (Mat. no. 4 692 440)

Mat. no. 4 837 600 Art. no. 69.7464.01

#### E 454 injector nozzle

- For trocar sleeves 10-15 mm, 8 x 150 mm
- Height-adjustment clip (Mat. no. 4 692 430)

Mat. no. 4 203 940 Art. no. 69.7454.01



#### E 471 spring clamp

For injector nozzle
 Ø 2.5 mm, for E 452
 Mat. no. 4 451 930
 Art. no. 69.7471.01

#### E 472 spring clamp

For injector nozzle
 Ø 4.0 mm, for E 351, E 453
 Mat. no. 4 451 960
 Art. no. 69.7472.01

### Cross-product accessories



**E 792**Silicone hose with LuerLock connector, male

Mat. no. 6 989 160 Art. no. 69.7792.01

#### E 790

Connector LuerLock adapters, female/ female

Mat. no. 6 390 710 Art. no. 69.7790.01

#### E 791 Connector

Connector

LuerLock female/male

Mat. no. 6 390 720 Art. no. 69.7790.01

### Adapters available individually:

- Pack of 4 x E 790 LuerLock female/female
- Pack of 4 x E 791
   LuerLock female/male
- Pack of 4 x E 792
   LuerLock male with silicone hose
   (160 mm)



#### A 6 cover net 1/2

- 215 x 445 mm
- Stainless-steel frame with polyamide mesh (particularly stable and durable)
- For example for E 142 mesh tray
- Suitable for ophthalmology

Mat. no. 7 217 650 Art. no. 69.5006.01



#### E 499 adapter

- Adapter for suction and irrigation cannulae with various connectors, and Verres cannulae
- Compatible with carriages with injector nozzles

Mat. no. 7 759 500 Art. no. 69.5499.01



#### E 142 insert 1/2

- DIN mesh tray
- 1 mm wire gauge
- Mesh size 5 mm
- 5 mm all-round frame
- 2 hinged handles
- Max. load 10 kg
- H 45/55, W 255, D 480 mm

Illustration shows E 142 with E 476 receptacles and E 477 stops

Mat. no. 3 830 300 Art. no. 69.5142.01



#### E 476 adapters

- For use in mesh trays with 5 mm mesh (e.g. E 142)
- 50 per bag
- For instruments with a diameter of 4 to 8 mm

Mat. no. 5 692 400 Art. no. 69.7476.01



#### E 479 adapters

- For use in mesh trays with 5 mm mesh (e.g. E 142)
- 50 per bag
- For instruments with a diameter of up to 4 mm

Mat. no. 6 553 750 Art. no. 69.7479.01



#### E 477 stoppers

- For use in mesh trays with 5 mm mesh (e.g. E 142)
- 20 per bag

Mat. no. 5 692 390 Art. no. 69.7477.01

### Inserts for GYN and ENT instruments



E 416 insert 1/4

- For 6 one or two-piece speculae
- 7 holders, spacing 40 mm
- For upper or lower basket
- H 157, W 178, D 279 mm

Mat. no. 3 888 650 Art. no. 69.5416.01



#### E 417/1 insert 2/5

- For approx. 30 ear and nose speculae and Brünings tongue spatulas
- 280 compartments, approx. 13 x 13 mm
- Mesh size on base 1.7 mm
- For upper or lower basket
- H 63, W 173, D 445 mm

Mat. no. 9 033 260 Art. no. 69.5417.02



#### E 803 insert 2/5

- For ear and nose speculae
- 160 compartments, approx. 13 x 13 mm
- Mesh size on base 1.7 mm
- For upper or lower basket
- H 63, W 165, D 317 mm

Mat. no. 4 637 750 Art. no. 69.5803.01



#### E 373 insert 1/6

- For ENT instruments (e.g. ear funnels)
- Wire mesh, mesh size: Base 3 mm, sides 1.7 mm, lid 3 mm
- 28 upright supports
- 2 hinged handles
- For upper or lower basket
- H 55, W 150, D 225 mm

Mat. no. 3 809 790 Art. no. 69.5373.01



#### E 374 insert 2/5

- For ENT instruments such as inhalation tubes, etc.
- 24 compartments, approx. 45 x 45 mm
- 27 compartments, approx. 12 x 12 mm
- Mesh size on base 1.7 mm
- For upper or lower basket
- H 63, W 173, D 445 mm

Mat. no. 3 809 800 Art. no. 69.5374.01

### Mesh trays/baskets



#### E 146 insert 1/6

- Mesh size on base 3 mm
- Mesh size on sides 1.7 mm
- Mesh size on lid 8 mm
- 2 hinged handles
- For upper or lower basket
- H 55, W 150, D 225 mm

Mat. no. 3 346 680 Art. no. 69.5146.01

#### E 363 insert 1/6

(not illustrated)

- Mesh size 1 mm, with lid
- For upper or lower basket
- H 55, W 150, D 225 mm

Mat. no. 3 149 790 Art. no. 69.5363.01



#### E 373 insert 1/6

- For ENT instruments (e.g. ear funnels)
- Wire mesh, mesh size:
   Base 3 mm, sides 1.7 mm,
   lid 3 mm
- 28 upright supports
- 2 hinged handles
- For upper or lower basket
- H 55, W 150, D 225 mm

Mat. no. 3 809 790 Art. no. 69.5373.01



#### E 441/1 insert 1/4

- For micro-instruments
- Mesh size on base 1.7 mm
- Solid sides, stackable
- Internal divisions with 6 adjustable supports provide the ideal storage for instruments.
- For upper or lower basket
- H 60, W 183, D 284 mm

Mat. no. 5 988 240 Art. no. 69.5441.02



#### E 131/1 insert 1/2

- For 5 mesh trays/kidney dishes
- 6 holders, H 160 mm, spacing 80 mm
- For lower basket
- H 168, W 180, D 495 mm

Mat. no. 5 771 030 Art. no. 69.5131.02



E 473/1 mesh insert with lid

- Mesh basket with lid for small parts
- Hooks in place
- H 85, W 60, D 60 mm

Mat. no. 5 341 570 Art. no. 69.5473.02



#### E 451 insert 1/6

- Mesh tray with lid for sundry small items
- Wire gauge:1 mm base0.8 mm sides,1 mm lidMesh gauge:3 mm base
- 1.7 mm sides
- 3 mm lid
- H 55, W 150, D 225 mm

Mat. no. 4 203 910 Art. no. 69.5451.01

### Cover nets, inserts, lids



E 379 insert 1/2

- For various utensils
- 0.8 mm wire gauge
- Mesh size 1.7 mm
- 5 mm all-round frame
- 2 handles
- For upper or lower basket
- H 80/110, W 180, D 445 mm

Mat. no. 3 203 500 Art. no. 69.5379.01



E 378 insert 1/1

- For various utensils
- 0.8 mm wire gauge
- Mesh size 1.7 mm
- 5 mm all-round frame
- 2 handles
- For lower basket
- H 80 (110), W 460, D 460 mm

Mat. no. 3 203 530 Art. no. 69.5378.01



E 143 insert 1/4

- Mesh tray
- 1 mm wire gauge
- Mesh size 5 mm
- 5 mm all-round frame
- 2 hinged handles
- Max. load 5 kg
- For upper or lower basket
- H 45 (55), W 255, D 230 mm

Mat. no. 3 830 310 Art. no. 69.5143.01



E 142 insert 1/2

- DIN mesh tray
- 1 mm wire gauge
- Mesh size 5 mm
- 5 mm all-round frame
- 2 hinged handles
- Max. load 10 kg
- H 45/55, W 255, D 480 mm

Mat. no. 3 830 300 Art. no. 69.5142.01



**UTS Utensil box** 

- For sundry small items, with lid
- 3 compartments, compartment size
   115 x 100 mm
- H 93, W 102, D 350 mm

Mat. no. 3 810 260 Art. no. 69.7410.01



E 430/1 insert 1/3 mesh tray

- Wire mesh, mesh size 5 mm
- H 40, W 150, D 445 mm

Mat. no. 6 643 200 Art. no. 69.7430.03



E 468 1/4 mesh insert with lid and compartments

- For various utensils
- Made from welded crosswires/stainless steel with lid and compartments
- Mesh size 5 x 5 x 1 mm
- H 70 (76), W 250, D 170 mm

Mat. no. 5 043 640 Art. no. 69.5468.01



#### A 18 1/2 lid

- For insert E 484
- For insert E 142 as 1/1 cover
- 1 mm wire mesh5 mm mesh gauge4 mm all-round frame
- H 20, W 430, D 230 mm

Mat. no. 5 631 080 Art. no. 69.7450.06

#### A 16 1/2 lid (not illustrated)

- For insert E 142
- 230 x 215 mm
- Stainless steel

Mat. no. 3 981 990 Art. no. 69.7450.04

### Inserts, mesh trays



#### E 908/1 insert

- For separable minimally invasive instruments or fibre optics
- Mesh width 8 x 1 mm (base), sides closed
- Internal dividers can be repositioned with 4 supports
- Hooks to slot into E 906/1
- H 36, W 130, D 460 mm

Mat. no. 5 704 010 Art. no. 69.5908.02



#### E 907/1 insert/mesh tray

- For sundry small items, with lid
- Mesh size 3 x 1 mm
- Hooks to slot into E 905/1 and E 903/1
- H 48, W 129, D 170 mm

Mat. no. 6 943 820 Art. no. 69.5907.02



#### E 460 insert 1/4

- For rigid fibre optics of varying lengths
- Mesh size, base8 x 1 mmSides/lid 7 x 7 x 3 mm
- With 3 holders for 2 rigid fibre optics of varying lengths
- H 53, W 100, D 430 mm

Mat. no. 4 861 390 Art. no. 69.5460.01



#### E 457 insert 1/2

- For separable MI instruments (e.g. 12 handles and inner shafts)
- Mesh base, 3 mm mesh size, solid sides
- Welded holders for 8-12 handles, variable arrangement, with 3 holders for storing and securing 8-12 shafts from disassembled MI instruments
- H 62, W 192, D 490 mm

Mat. no. 4 618 030 Art. no. 69.5457.02



#### E 473/1 insert/filter

- Mesh basket with lid for small parts
- For hooking onto mesh trays
- H 85, W 60, D 60 mm

Mat. no. 5 341 570 Art. no. 69.5473.02



#### E 444 insert/spiral rack

- For fibre optics and irrigation tubes
- Fibre optics and tubes wound spirally around central rack
- H 168 mm, with folding clip, 214 mm
- •Ø 140 mm

Mat. no. 4 637 700 Art. no. 69.5444.01



E 362 blanking screw

 M 8 x 1 thread, to close connectors on load carriers

Mat. no. 3 809 630 Art. no. 69.7362.01

# Transport trolleys



#### MF/3 for G 7823/G 7824/G 7825/G 7826

- Trolley to simplify handling of mobile injector units
- Footswitch-operated height-adjustment mechanism
- 4 wheels, of which 2 are lockable
- H 1182, W 660, D 807 mm, +/- 100 mm
- Docking height 751 mm, +/- 100 mm

Mat. no. 6 392 900 Art. no. 69.2001.07



#### MF 27/28-1 for PG 8527/PG 8528

- 4 wheels, of which 2 are lockable
- Both ends can dock onto machine, tables, sluices or conveyors.
- Docking height 850, -100, + 150 mm
- H 1050, W 740, D 930 mm
- With removable drip tray H 70, W 603, D 866 mm

Mat. no. 7 397 640 Art. no. 69.2001.11





### Washer-disinfectors G 7823 und G 7824

The G 7823 and G 7824 washer-disinfectors were specifically designed to meet the needs of smaller and mediumsized hospitals and outpatient surgeries. With a width of only 900 mm, these machines are the ideal proposition in cases where space is at a premium. Depending on the installation site and hygiene requirements, machines are available as single-door front-loading models (G 7823) or as two-door barrier machines (G 7824) for installation in a diaphragm wall.

### Flexible solutions for the central and decentral reprocessing of instruments

Miele's G 7823 and G 7824 washerdisinfectors offer a wide range of installation options, tailored individually to the needs of hospitals and outpatient surgeries. This allows units to be installed both **decentrally** in specialist departments or centrally in **CSSD** for coping with larger volumes of instruments.

| Performance comparison<br>Washer-disinfectors | Unit width/depth | Doors                 | Cabinet Useable dimensions H/W/D Useable volume | Approx.<br>capacity per cycle                |
|---|------------------|-----------------------|---|--|
| PG 8527 / PG 8528                             | 1150/870 mm      | Vertical-rise door(s) | 675/650/800 mm<br>351 I                         | 18 DIN mesh trays<br>7 AN sets<br>3 MI sets  |
| G 7823 / G 7824                               | 900/770 mm       | Bottom-hinged door(s) | 510/530/620 mm<br>168 l                         | 8 DIN mesh trays<br>3 AN sets:<br>2 MI sets  |
| G 7825 / G 7826                               | 900/750 mm       | Bottom-hinged door(s) | 683/541/610 mm<br>225 l                         | 10 DIN mesh trays<br>4 AN sets:<br>2 MI sets |

### Machine features and specifications



#### **Versions**

- G 7823: Front-loading unit with single bottom-hinged door
- G 7824: Two-door barrier washerdisinfector for installation in a diaphragm wall

#### Design

- Stand-alone or side-by-side installation
- Width 900 mm
- Modular approach with customised features to meet individual requirements
- Single-chamber system for washing, disinfection and drying
- Service-friendly design
- Heater elements outside wash chamber
- Low heat and sound emissions thanks to double insulation

#### Capacity per cycle

or 2 MI sets

8 DIN mesh trays
or 3 DIN containers
150 x 300 x 600 mm
or 1 DIN container
300 x 300 x 600 mm
and 1 DIN container 150 x 300 x 600 mm
or 3 AN sets

#### Cleaning technology

- Hygienic freshwater system with fresh water intake in each programme stage
- Cleaning, disinfection and drying in a closed, single-cabinet system
- 2 spray arms in cabinet for thorough cleaning of instrument surfaces
- Spray arms with high water impact force
- Full water jet access, ensuring optimum results
- Thorough cleaning of lumens thanks to injector system
- Direct docking of load carriers to water circuit



#### Standard features

- 2 powerful circulation pumps
- Triple filtration with large surface filter, coarse filter and micro-fine filter
- Filters in inlet hoses
- Flowmeter to monitor water intake quantities
- 1 dump valve

#### **Dispenser systems**

- 2 dispenser pumps for liquid detergent and acidic agents
- Extendible (further 3 internal dispenser pumps retrofittable)

#### Water softeners

 High-capacity water softener (integratable)

#### Drainage

• 2 drain pumps

#### Controls

- Freely programmable controls PROFITRONIC
- 64 programme slots
- 13 standard washing and disinfection programmes
- 8 service programmes
- 43 vacant programme slots
- User interface with local-language display
- Display of programme selection and programming dialogs, countdown time, faults and operating hours.
- Compilation of new programmes using machine controls or using PC/laptop via optical interface

#### Interfaces

- RS 232 serial interface for process documentation and remote service
- Optical interface for service and maintenance

#### **Printer**

 Integrated 6-pin printer to record key process parameters

#### Safety features

- Electric door lock
- Programme recontinuation in event of power outage
- Peak-load negotiation
- Optical and acoustic signal at end of programme
- 2 separate sensors for temperature monitoring and control
- Port for simple positioning of sensors in the wash cabinet for process validation
- Sensors in cabinet and a magnetic strip on load carriers for automatic assignment of programmes to load carriers

# Modular machine concept Optional extras





#### SBWR/3

#### Roller plinth/floor tray for G 7823

- Roller plinth/floor tray
- Castors allow machine to be pulled forward for servicing
- H 100, W 900, D 768 mm

#### TA-S1

#### Drying unit (optional)

- Hot-air drying unit (electrically heated)
- Low-maintenance three-phase motor with side-channel compressor
- Temperature range 60–115°C
- Time range 1-240 mins.
- 1 x pre-filter EU 4, Filter rating > 95% Filter life 200 h
- 2 particulate filters H 13, Filter rating > 99.992%, Filter life 500 h



#### SBW/4

#### Plinth/floor tray

- Frame with built-in floor tray
- Cut-outs for supplies and utilities
- A continuous plinth/floor tray can be provided on site when multiple machines are installed in a row.
- For G 7824
- H 100, W 900, D 768 mm

#### DK 23/24 steam condenser

- Steam condenser (heat exchanger)
- No water consumption when connected to on-site cooling water circuit, max. water pressure 8 bar

# Machine versions, components, accessories

#### **MAV 23/24**

#### **Top-box panelling**

- Stainless-steel panelling with lockable service hatches to conceal steam condenser
- MAV without cover panel
- Panelling between top box and ceiling must be provided on site
- H 512, W 900, D 768 mm

#### Cover plate

• Top cover for MAV 23/24

#### DE 23 lid

 Lid for G 7823 when installed without steam condenser (not suitable for G 7824)

#### **Retrofit kits**

- OxiVario (OXI-NRS 1-G 7823/24)
- OrthoVario (ORTHO)

Retrofit kits and on-site installation by Miele Service

#### Note

An emergency door release must be fitted if the G 7824 is installed without top panelling (MAV). The required parts are available from Spares (cf. Conversion Instructions, Mat. no. 6476390).

| G 7823                                    | Features  | Mat. no. | Art. no.   |
|---|---|----------|------------|
| Electric                                  | AE PT EL AP                                     | 6174400  | 62.7823.24 |
| Steam                                     | AE PT D AP                                      | 6174420  | 62.7823.14 |
| Steam/Electric                            | AE PT D/EL AP                                   | 6174430  | 62.7823.34 |
| G 7824                                    |   |          |            |
| Electric                                  | AE PT EL AV                                     | 6174350  | 62.7824.21 |
| Electric                                  | AE PT EL AP                                     | 6174360  | 62.7824.24 |
| Steam                                     | AE PT D AV                                      | 6174370  | 62.7824.11 |
| Steam/Electric                            | AE PT D/EL AV                                   | 6174380  | 62.7824.31 |
| Accessories                               |   |          |            |
| WES-GE                                    | Water softeners                                 | 6174450  | 69.2020.01 |
| DOS 10/30                                 | Dispenser - pre-installed                       | 6224670  | 69.2250.04 |
| DOS 60/30                                 | Dispenser - pre-installed                       | 6224680  | 69.2250.05 |
| PRT/1                                     | Printer - pre-installed                         | 6224720  | 69.2211.03 |
| OxiVario module                           | OxiVario – pre-installed                        | 6444140  | 69.2250.06 |
| OrthoVario module (incl. OxiVario module) | OrthoVario - pre-installed                      | 6946150  | 69.2250.08 |
| ML  | 5 magnets*                                      | 5267440  | 69.7457.02 |
| SBWR/3                                    | Roller plinth/floor tray                        | 6174480  | 69.3710.06 |
| SBW/4                                     | Plinth/floor tray                               | 6174490  | 69.3710.07 |
| TA-S1                                     | Integrated drying unit                          | 6174440  | 69.2010.03 |
| DK 23/24                                  | Steam condenser                                 | 6600610  | 69.2300.06 |
| MAV 23/24                                 | Top-box panelling                               | 6805170  | 69.2100.09 |
| Cover plate                               | For MAV 23/24                                   | 6805180  | 69.2100.10 |
| DE 23                                     | Lid for G 7823                                  | 6174470  | 69.2161.01 |
| OXI-NRS 1-G 7823/24                       | OxiVario retrofit kit, including installation   | 6444690  | 69.2310.05 |
| ORTHO (excl. OxiVario)**                  | OrthoVario retrofit kit, including installation | 6930680  |            |

Explanation of symbols cf. Page 5

<sup>\*</sup> for automatic load carrier recognition. Not included with load carriers for G 7823/G 7824.

<sup>\*\*</sup> presumes installation of OxiVario

# Technical data

| Front-loading unit with bottom-hinged door Barrier model with bottom-hinged doors Single/multiple installations Freshwater system, max. temperature 93°C Direct load carrier docking for cleaning and drying of lumened instruments 2 circulation pumps [Qmax. l/min] Controls/Programmes | •<br>-<br>• | -        |
|---|-------------|----------|
| Single/multiple installations reshwater system, max. temperature 93°C Direct load carrier docking for cleaning and drying of lumened instruments circulation pumps [Qmax. I/min]  | •           |          |
| reshwater system, max. temperature 93°C<br>Direct load carrier docking for cleaning and drying of lumened instruments<br>circulation pumps [Qmax. I/min]  |             | •        |
| Direct load carrier docking for cleaning and drying of lumened instruments 2 circulation pumps [Qmax. I/min]  | •           | •        |
| 2 circulation pumps [Qmax. I/min]   |             | •        |
|   | •           | •        |
| Controls/Programmes   | 300/500*    | 300/500* |
|   |             |          |
| Profitronic, freely programmable  | •           | •        |
| 64 programme slots  | •           | •        |
| Electric door lock  | •           | •        |
| Peak-load negotiation   | •           | •        |
| Serial interface for process documentation  | •           | •        |
| Magnetic strip for automatic load carrier recognition   | •           | •        |
| Remote service enabled  | •           | •        |
| Water connections   |             |          |
| 1 x Cold water, 2-10 bar flow pressure (200-1000 kPa) 4°dH  | •           | •        |
| 1 x Hot water, 2–10 bar flow pressure (200–1000 kPa) 4°dH   | •           | •        |
| 1 x demineralised water, 2–10 bar flow pressure (200–1000 kPa)  | •           | •        |
| 3 inlet hoses ½" with ¾" threaded union   | •           | •        |
| Drain valve DN 50, odour trap to be fitted on site  | •           | •        |
| 2 drain pumps, odour trap to be provided on site, drain hose DN 22  | 0           | 0        |
| Electrical connection: Electric heating   |             |          |
| 3 N AC 400 V 50 Hz  | •           | •        |
| Heating [kW]  | 9.0         | 9.0      |
| Circulation pump [kW]   | 0.4/0.8     | 0.4/0.8  |
| Total rated load w/o TA [kW]  | 10.2        | 10.2     |
| Total rated load with electric TA [kW]  | 10.2        | 10.2     |
| Fuse rating [A]   | 3 x 16      | 3 x 16   |
| Electrical connection: Steam heating  | 0 X 10      | 0 X 10   |
| 3 N AC 400 V 50 Hz  | •           | •        |
| Circulation pump [kW]   | 0.4/0.8     | 0.4/0.8  |
| Total rated load w/o TA [kW]  | 1.5         | 1.5      |
| Total rated load with electric TA [kW]  | 4.5         | 4.5      |
| Fuse rating [A]   | 3 x 16      | 3 x 16   |
| Steam connection G ½" (DN 15)   | •           | •        |
| Operating pressure 250–1000 kPa   | •           | •        |
|   | •           | •        |
| Compressed air connection 600–1200 kPa  |             | •        |
| Electrical connection: Steam/Electric (convertible)   | •           | •        |
| 3 N AC 400 V 50 Hz  |             |          |
| Heating [kW]  | 9.0         | 9.0      |
| Circulation pump [kW]   | 0.4/0.8     | 0.4/0.8  |
| Total rated load with electric TA [kW]  | 10.2        | 10.2     |
| Fuse rating [A]   | 3 x 16      | 3 x 16   |
| Steam connection G ½" (DN 15)   | •           | •        |
| Operating pressure 250–1000 kPa   | •           | •        |
| Compressed air connection 600–1200 kPa  | •           | •        |

# Technical data

| ### Spenser systems  ### A DOS 1070 dispenser pump for liquid acidic agents  ### A DOS 8070 dispenser pump for liquid acidic agents  ### A 15 supply caristers  ### A 15 supply caristers  ### A 20 Supply cariste | Washer-disinfectors | G 7823 (1 door) | G 7824 (2 door) |
|--|---------------------|-----------------|-----------------|
| x DOS 10/30 dispenser pump for liquid acidic agents         •         •           x DOS 60/30 dispenser pump for liquid detergent         •         •           x 5 I supply canisters         •         •           pace for 2 x 5 I supply canisters         •         •           pace for 2 x 5 I supply canisters         •         •           urther connection options, internal pumps (optional)         v         v           x DOS 10/30 dispenser pump for surfactant and neutralising agent         •         •           x DOS 60/30 dispenser pump for detergent or disinfectant         •         •           an [kW]         0.9         0.9           eater bank [kW]         3.6         3.6           otal rated load [kW]         4.5         4.5           int throughput [m³/h]         125         125           int perature selection in 1° increments         60-115         60-115           int e selection in 1-minute increments         1-240         1-240           int perature selection in 1. filter rating > 95%, filter life 200 h         •         •           intensions/Weight         v         v           intensions Wight         v         v           intensions Hincl. plinth/tray and MAV 23/24 (top panelling installation kit) [mm]         1928  |                     | G 7020 (1 0001) | G 1027 (2 0001) |
| x DOS 60/30 dispenser pump for liquid detergent         •         •           x 5 I supply canisters         •         •           pace for 2 x 5 I supply canisters         •         •           wither connection options, internal pumps (optional)         .         .           x DOS 10/30 dispenser pump for surfactant and neutralising agent         •         •           x DOS 60/30 dispenser pump for detergent or disinfectant         •         •           an [kW]         0.9         0.9           acater bank [kW]         3.6         3.6           obat rated load [kW]         4.5         4.5           in throughput [m³/h]         125         125           interesture selection in 1° increments         60–115         60–115           interesture selection in 1° increments         1–240         1–240           interesture selection in 1° increments         9.         1           interesture selection in 1° increments         9.         9.           interesture selection in 1° increments         9.         9.   |                     | •               | •               |
| x 5 I supply canisters         •         •           pace for 2 x 5 I supply canisters         •         •           urther connection options, internal pumps (optional)         o         o           x DOS 10/30 dispenser pump for surfactant and neutralising agent         o         o           x DOS 60/30 dispenser pump for detergent or disinfectant         o         o           rying unit/side-channel compressor         u         0           an [kW]         0.9         0.9           otal rated load [kW]         3.6         3.6           otal rated load [kW]         4.5         4.5           in throughput [m³/h]         125         125           eater bank [ker pating selection in 1° increments         60–115         60–115           xi pre-flitters EU 4, filter rating > 95.96, filter life 200 h         •         •           xi pre-flitters EU 4, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h         •         •           xi pre-flitters EU 4, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h         •         •           xi pre-flitters EU 4, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h         •         •           xi pre-flitters EU 4, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h         •         •           xi pre-flitters EU 4, filter r  |                     | •               | •               |
| Pace for 2 x 5   supply canisters   Pace for 2 x 5   supply cani   |                     | •               | •               |
| arther connection options, internal pumps (optional)         Control 10/30 dispenser pump for surfactant and neutralising agent         o         o           x DOS 60/30 dispenser pump for detergent or disinfectant         o         o         o           x DOS 60/30 dispenser pump for detergent or disinfectant         o         o         o           x DOS 60/30 dispenser pump for detergent or disinfectant         o         o         o           ar (kW)         0.9         0.9         0.9           deter bank [kW]         3.6         3.6         3.6           obal rated load [kW]         4.5         4.5         4.5           introluphut [m²/n]         125         125         125           sumperature selection in 1° increments         60–115         60–115         60–115           sumperature selection in 1-minute increments         1–240         1–240         1–240           x pre-filters EU 4, filter rating > 95%, filter life 200 h         •         •         •           sumperature selection in 1-minute increments         1–240         •         •           sternal dimensions Meight         90%         •         •         •           sternal dimensions William         90%         90%         8         8         8         8         8   |                     | •               | •               |
| x DOS 10/30 dispenser pump for surfactant and neutralising agent         o         o           x DOS 60/30 dispenser pump for detergent or disinfectant         o         o           x DOS 60/30 dispenser pump for detergent or disinfectant         o         o           x DOS 60/30 dispenser pump for detergent or disinfectant         o         o           x POS 10/30 dispenser pump for detergent or disinfectant         o         o           x POS 60/30 dispenser pump for detergent or disinfectant         o         o           x POS 10/30 dispenser pump for detergent or disinfectant         o         o           x POS 10/30 dispenser pump for detergent or disinfectant         o         o           x POS 10/30 dispenser pump for detergent or disinfectant         o         0           x POS 10/30 dispenser pump for detergent or disinfectant         0         0         0           acter plant dispenser pump for detergent or disinfectant         4.5         4.5         4.5           ster tribugh determined of the pank [kW]         4.5         4.5         4.5         4.5         4.5         4.5         4.5         4.5         4.5         4.5         4.5         4.5         4.5         4.2         4.0         4.2         4.0         4.2         4.0         4.0         4.0         4.0         4.0  |                     |                 |                 |
| x DOS 60/30 dispenser pump for detergent or disinfectant         o         o           rying unit/side-channel compressor         w         v           an [kW]         0.9         0.9           eater bank [kW]         3.6         3.6           otal rated load [kW]         4.5         4.5           otal rated load [kW]         60-115         60-115           in throughput [m³/h]         125         125           sime selection in 1° increments         60-115         60-115           time selection in 1-minute increments         1-240         1-240           x pre-filters EU 4, filter rating > 95%, filter life 200 h         0         0           x pre-filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h         0         0           sternal dimensions Wingth         900/768         1928         1928           sternal dimensions H incl. plinth/tray and MAV 23/24 (top panelling installation kit) [mm]         1928         1928         1928           sternal dimensions Wingth         900/768         900/768         900/768         900/768         900/768         900/768         900/768         900/768         900/768         900/768         900/768         900/768         900/768         900/768         900/768         900/768         900/   |                     | 0               | 0               |
| rying unit/side-channel compressor an [kW] 0.9 0.9 aeater bank [kW] 3.6 3.6 bal rated load [kW] 4.5 4.5 ir throughput [m³/h] 125 125 imperature selection in 1° increments 60–115 ime selection in 1-minute increments 1–240 1–240 x pre-filters EU 4, filter rating > 95%, filter life 200 h particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h intensions/Weight  sternal dimensions H incl. plinth/tray and MAV 23/24 (top panelling installation kit) [mm] 1928 1928 xternal dimensions W/D [mm] 900/768 seable cabinet dimensions H/W/D [mm] 700/530/620 510/530/620 verall cabinet dimensions H/W/D [mm] 865 865 leight [kg] 350 370 seaing tainless steel (AE) • • • • • • • • • • • • • • • • • • •   |                     |                 |                 |
| an [kW]       0.9       0.9         eater bank [kW]       3.6       3.6         batal rated load [kW]       4.5       4.5         in throughput [m³/h]       125       125         emperature selection in 1° increments       60–115       60–115         ime selection in 1-minute increments       1–240       1–240         x pre-filters EU 4, filter rating > 95%, filter life 200 h       •       •         particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h       •       •         imensions/Weight       4       5       900/768         xternal dimensions H incl. plinth/tray and MAV 23/24 (top panelling installation kit) [mm]       1928       1928         xternal dimensions WD [mm]       900/768       900/768         seable cabinet dimensions H/W/D [mm]       700/530/620       510/530/620         verall cabinet dimensions H/W/D [mm]       700/530/660       700/530/600         ocking height above floor (with plinth) [mm]       865       865         verall cabinet selection in 1-minute increments       900/768       865         verall cabinet dimensions H/W/D [mm]       700/530/660       700/530/660         ocking height above floor (with plinth) [mm]       865       865         verall cabinet dimensions H/W/D   |                     |                 |                 |
| eater bank [kW]       3.6       3.6         bala rated load [kW]       4.5       4.5         in throughput [m³/h]       125       125         emperature selection in 1° increments       60–115       60–115         ime selection in 1-minute increments       1–240       1–240         x pre-filters EU 4, filter rating > 95%, filter life 200 h       •       •         particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h       •       •         imensions/Weight       xternal dimensions H incl. plinth/tray and MAV 23/24 (top panelling installation kit) [mm]       1928       1928         xternal dimensions W/D [mm]       900/768       900/768         seable cabinet dimensions H/W/D [mm]       510/530/620       510/530/620         verall cabinet dimensions H/W/D [mm]       865       865         verall cabinet dimensions H/W/D [mm]       865       865         verall peight above floor (with plinth) [mm]       865       865         verall cabinet selection in 1-minute increments       1-240       1-240         xternal dimensions H/W/D [mm]       510/530/620       700/530/620         verall cabinet dimensions H/W/D [mm]       865       865         verall cabinet dimensions H/W/D [mm]       865       865  |                     | 0.9             | 0.9             |
| Stal rated load [kW]   4.5     |                     |                 |                 |
| ir throughput [m³/h]       125       125         emperature selection in 1° increments       60–115       60–115         ime selection in 1-minute increments       1–240       1–240         x pre-filters EU 4, filter rating > 95%, filter life 200 h       •       •         particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h       •       •         imensions/Weight       verall dimensions H incl. plinth/tray and MAV 23/24 (top panelling installation kit) [mm]       1928       1928         external dimensions W/D [mm]       900/768       900/768         seable cabinet dimensions H/W/D [mm]       510/530/620       510/530/620         everall cabinet dimensions H/W/D [mm]       700/530/660       700/530/660         ocking height above floor (with plinth) [mm]       865       865         /eight [kg]       350       370         asing         tainless steel (AE)       •       •         extertificates         DE, VDE-EMC, MDD CE 0366       •       •   |                     |                 |                 |
| emperature selection in 1° increments       60–115       60–115         ime selection in 1-minute increments       1–240       1–240         ix pre-filters EU 4, filter rating > 95%, filter life 200 h       •       •         particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h       •       •         imensions/Weight       verall dimensions H incl. plinth/tray and MAV 23/24 (top panelling installation kit) [mm]       1928       1928         ixternal dimensions W/D [mm]       900/768       900/768         seable cabinet dimensions H/W/D [mm]       510/530/620       510/530/620         verall cabinet dimensions H/W/D [mm]       700/530/660       700/530/660         ocking height above floor (with plinth) [mm]       865       865         deight [kg]       350       370         asing         tainless steel (AE)       •       •         ost certificates         DE, VDE-EMC, MDD CE 0366       •       •  |                     |                 |                 |
| ime selection in 1-minute increments  x pre-filters EU 4, filter rating > 95%, filter life 200 h particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h particulate filters/HEPA H 2, so  |                     |                 |                 |
| x pre-filters EU 4, filter rating > 95%, filter life 200 h particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h imensions/Weight  xternal dimensions H incl. plinth/tray and MAV 23/24 (top panelling installation kit) [mm] 1928 1928  xternal dimensions W/D [mm] 900/768 900/768 seable cabinet dimensions H/W/D [mm] 510/530/620 510/530/620 verall cabinet dimensions H/W/D [mm] 700/530/660 700/530/660 ocking height above floor (with plinth) [mm] 865 865 leight [kg] 350 370  asing  tainless steel (AE) • • • • • • • • • • • • • • • • • • •  |                     |                 |                 |
| particulate filters/HEPA H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h  imensions/Weight  xternal dimensions H incl. plinth/tray and MAV 23/24 (top panelling installation kit) [mm] 1928 1928  xternal dimensions W/D [mm] 900/768 900/768  seable cabinet dimensions H/W/D [mm] 510/530/620 510/530/620  verall cabinet dimensions H/W/D [mm] 700/530/660 700/530/660  ocking height above floor (with plinth) [mm] 865 865  /eight [kg] 350 370  asing  tainless steel (AE) • • • • • • • • • • • • • • • • • • •  |                     |                 |                 |
| imensions/Weight         xternal dimensions H incl. plinth/tray and MAV 23/24 (top panelling installation kit) [mm]         1928         1928           xternal dimensions W/D [mm]         900/768         900/768           seable cabinet dimensions H/W/D [mm]         510/530/620         510/530/620           verall cabinet dimensions H/W/D [mm]         700/530/660         700/530/660           ocking height above floor (with plinth) [mm]         865         865           /eight [kg]         350         370           asing         *         *           tainless steel (AE)         *         *           est certificates         *         *           DE, VDE-EMC, MDD CE 0366         *         *   |                     |                 |                 |
| Section   Sect   |                     | •               | •               |
| external dimensions W/D [mm]       900/768       900/768         seable cabinet dimensions H/W/D [mm]       510/530/620       510/530/620         everall cabinet dimensions H/W/D [mm]       700/530/660       700/530/660         ocking height above floor (with plinth) [mm]       865       865         /eight [kg]       350       370         asing       ***       ***         tainless steel (AE)       ***       ***         est certificates       ***       ***         DE, VDE-EMC, MDD CE 0366       ***       **  |                     | 1029            | 1020            |
| seable cabinet dimensions H/W/D [mm]       510/530/620       510/530/620         verall cabinet dimensions H/W/D [mm]       700/530/660       700/530/660         ocking height above floor (with plinth) [mm]       865       865         /eight [kg]       350       370         asing         tainless steel (AE)       •       •         est certificates         DE, VDE-EMC, MDD CE 0366       •       •   |                     |                 |                 |
| verall cabinet dimensions H/W/D [mm]       700/530/660       700/530/660         ocking height above floor (with plinth) [mm]       865       865         /eight [kg]       350       370         asing       **       **         tainless steel (AE)       **       **         est certificates       **       *         DE, VDE-EMC, MDD CE 0366       *       *   | • •                 |                 |                 |
| ocking height above floor (with plinth) [mm]       865       865         /eight [kg]       350       370         asing       *       *         tainless steel (AE)       *       *         est certificates       *       *         DE, VDE-EMC, MDD CE 0366       *       *   |                     |                 |                 |
| /eight [kg]     350     370       asing     ***     ***       tainless steel (AE)     ***     ***       est certificates     ***     ***       DE, VDE-EMC, MDD CE 0366     **     **  |                     |                 |                 |
| asing         •         •           tainless steel (AE)         •         •           est certificates         •         •           DE, VDE-EMC, MDD CE 0366         •         •  |                     |                 |                 |
| tainless steel (AE)  est certificates  DE, VDE-EMC, MDD CE 0366  • • •   |                     | 350             | 370             |
| est certificates DE, VDE-EMC, MDD CE 0366  • • •   |                     | _               | _               |
| DE, VDE-EMC, MDD CE 0366 • •   |                     | •               | •               |
|  |                     |                 |                 |
| • • • • • • • • • • • • • • • • • • •  |                     | •               | •               |
|  |                     |                 |                 |
|  |                     |                 |                 |
|  |                     |                 |                 |
|  |                     |                 |                 |
|  |                     |                 |                 |
|  |                     |                 |                 |
|  |                     |                 |                 |
|  |                     |                 |                 |
|  |                     |                 |                 |
|  |                     |                 |                 |
|  |                     |                 |                 |
|  |                     |                 |                 |
|  |                     |                 |                 |
|  |                     |                 |                 |
|  |                     |                 |                 |
|  |                     |                 |                 |
| = Standard, o = Optional   |                     |                 |                 |

### Load carrier with 1 - 3 levels



#### E 550 load carrier

- Load carrier with 1 level
- For use with various inserts
- Magnetic strip for automatic load carrier recognition (excl. ML magnets)
- H 250, W 530, D 645 mm
- Loading dimensions:
   H 481, W 482, D 587 mm
   Mat. no. 6 415 970

Art. no. 69.5550.01

#### E 550 load carrier with E 484 insert

(no illustration, description on P. 32)

• For various utensils



#### E 550 with E 730 insert

- Load carrier and insert for 20 theatre shoes
- Magnetic strip for automatic load carrier recognition (excl. ML magnets)
- H 298, W 530, D 645 mm

#### E 730

(no illustration, description on P. 32) Mat. no. 5 271 020 Art. no. 69.5730.01



#### E 484 insert with E 487

- E 484 can be fitted with 4 x E 487 holders each for 4 theatre shoes Height 280 mm
- Dimensions E 487 H 280, W 464, D 10 mm

#### E 484

Mat. no. 5 109 190 Art. no. 69.5484.01

#### E 487

(Description on Page 32) Mat. no. 5 109 250 Art. no. 69.5487.01



### E 525/2 mobile injector unit with drying

- Mobile injector unit with 3 levels, extendible to 4 levels (with E 702 module)
- Vertical clearance (bottom to top): 80/80/219 mm
- Magnetic strip for automatic load carrier recognition (excl. ML magnets)
- H 298, W 530, D 645 mm

Mat. no. 7 460 260 Art. no. 69.5525.03



E 525/2 load carrier with E 702 module

- Load carrier with 4 levels
- Vertical clearance (bottom to top): 80/80/90/80 mm
- Magnetic strip for automatic load carrier recognition (excl. ML magnets)
- H 458, W 530, D 645 mm

Illustration shows sample load with 8 x E 142 DIN mesh trays



#### E 702 top module

- Additional level for E 701/2 and E 525/2 load carriers
- For 2 DIN mesh trays
- H 160, W 530, D 560 mm

Mat. no. 5 221 490 Art. no. 69.5702.01



• Max. load per mesh tray 10 kg

# Load carriers with 3–4 levels Load carrier for containers



#### E 524/2 load carrier

- Load carrier with 3 levels
- Vertical clearance per level 125 mm
- Magnetic strip for automatic load carrier recognition (excl. ML magnets)
- H 390, W 530, D 645 mm

Illustration shows sample load with 6 x E 142 DIN mesh trays Mat. no. 7 460 250 Art. no. 69.5524.03

#### E 750 insert frame (not illustrated)

Mat. no. 5 567 510 Art. no. 69.5750.01



E 703/1 top module for E 525/2 and E 701

- Injector module for 7 lumened instruments
- H 170, W 80, D 555 mm
   (cf. P. 31 for items included)

Mat. no. 7 025 350 Art. no. 69.5703.02



#### E 555/2 load carrier

- Load carrier with 2 levels
- Vertical clearance (bottom to top): 210/219 mm
- Lower level W 487, D 592 mm
- Upper level W 495, D 548 mm
- Magnetic strip for automatic load carrier recognition (excl. ML magnets)
- H 298, W 530, D 645 mm

Mat. no. 7 460 290

Art. no. 69.5555.03

Illustration shows sample load with E 135 containers for baby bottles and E 364 inserts for teats.

E 750 insert frame required for E 135 container for baby bottles.



E 527/2 load carrier

- Load carrier for 1 DIN container and lid 300 x 600 mm and 150 x 600 mm
- Magnetic strip for automatic load carrier recognition (excl. ML magnets)
- H 500, W 530, D 645 mm

Mat. no. 7 460 270 Art. no. 69.5527.03



E 936 insert for boom lamp handles

- For various lamp handles
- 15 receptacles
- L 500, W 245, H 90 mm

Mat. no. 7 759 480 Art. no. 69.5936.01



E 528/2 load carrier

- Load carrier for 3 small DIN containers and lids, 150 x 600 mm
- Magnetic strip for automatic load carrier recognition (excl. ML magnets)
- H 500, W 530, D 645 mm

Mat. no. 7 460 280 Art. no. 69.5528.03

### Load carriers for MI instruments





#### E 499 adapter

- Adapter for suction and irrigation cannulae with various connectors, and Verres cannulae
- Compatible with carriages with injector nozzles

Mat. no. 7 759 500 Art. no. 69.5499.01

### E 526/1 mobile injector unit with drying connector

- Load carrier with 2 levels for MI instruments and accessories
- Vertical clearance (bottom) 100 mm
- Max. length of MI instruments: 600 mm
- Connection for hot-air drying unit
- Magnetic strip for automatic load carrier recognition (excl. ML magnets)
- H 498, W 530, D 645 mm

Illustration shows E 526/1 with E 451 mesh tray for sundry small items, insert for fibre-optics, E 457 mesh tray for separable MI instruments.

#### Supplied as standard with:

- 3 x E 336 injector sleeves MIBO for pipettes/MI instruments, 11 mm
- 2 x E 362 blanking screws
- 15 x E 442 injector sleeves for MI instruments, 4–8 mm
- 5 x E 443 injector sleeves for MI instruments, 8–8.5 mm
- 1 x E 444 insert/spiral rack for fibre optic cables or tubing
- 1 x E 445, 12 caps for irrigation sleeves/ MI instruments, 6 mm
- 1 x E 446, 12 caps for irrigation sleeves/ MI instruments, 10 mm
- 3 x E 447 adapters, female for male LuerLock
- 6 x E 448 silicone hoses, 300 mm, 5 x 1.5 mm
- 5 x E 449 adapters, male, for female LuerLock

- 1 x E 451 insert 1/6 mesh tray, H 55, W 150, D 225 mm
- 3 x E 452 injector nozzles, 2.5 x 60 mm
- 8 x E 453 injector nozzles,
  4.0 x 110 mm with spring supports
- 6 x E 454 injector nozzles for trocar sleeves, 10-15 mm
- 4 x E 456 spring stays for MI instruments
- 3 x E 464 sleeves for E 454 injector nozzles
- 2 x E 472 clamp for injector nozzle,
  4.0 mm
- 2 m silicone hose, Ø 5 mm Part no. 4822830

Mat. no. 7 026 560 Art. no. 69.5526.02

# Load carrier for MI instruments and reprocessing dental transmission instruments



#### E 574 load carrier with drying connector

- Modular system for approx. 2 MI OP sets
- For lumened instruments in 3 modules with integrated nozzles/adapters (not included)
- Applications:
- E 903/1 module for MI instruments/ urology
- E 905/1 module for short MI instruments E 906/1 module for long MI instruments 3 x E 444 drums for lightguides and irrigation hoses
- 1 x E 460 insert for rigid fibre optics 1 x E 457 insert for separable MI instruments or
- 1 x E 142 DIN mesh tray
- Connection for hot-air drying unit
- Magnetic strip for automatic load carrier recognition (excl. ML magnets)
- H 500, W 530, D 645 mm

#### Supplied as standard with:

- 3 x E 447 female adapters, for male LuerLock
- 6 x E 362 blanking screws

Mat. no. 6 802 840 Art. no. 69.5574.01

### Load carriers for micro-instruments



Geuder AG recommends that their current range of instruments be reprocessed in a Miele washer-disinfector.



#### E 529 load carrier with drying connector

- For 4–5 instrument sets (cataracts)
- Load carrier with 3 levels for microinstruments and accessories:
- 2 levels for mesh trays
- 1 level with 16 injector nozzles (hose and male LuerLock adapter) or
   30 injector nozzles (15 male and 15 female LuerLock adapters)
- Connection for E 478 on feed pipe
- Vertical clearance (bottom to top): 124/114.5/160 mm
- Connection for hot-air drying unit
- Magnetic strip for automatic load carrier recognition (excl. ML magnets)
- H 471, W 530, D 645 mm

Illustration shows E 529 with 4 x E 142 DIN mesh trays, with E 476 receptacles and E 477 end stoppers for micro-instruments.

Mat. no. 6 415 950 Art. no. 69.5529.01



#### FP ED

- Stainless-steel filter insert for E 478/1
- Diameter 30 mm

Mat. no. 7 290 770 Art. no. 69.7512.02



#### E 478/1 adapter

For 4 narrow-lumened cannulae (Sautter cannulae)

Mat. no. 7 292 920 Art. no. 69.7478.02

### Load carriers for anaesthetic equipment



### E 515/1 mobile injector unit with drying connector

- For approx. 3 AN sets or 9 breathing tubes, max. length 1.5 m, and various accessories
- Tubes wound around spiral drum
- E 430 mesh tray
- Connection for hot-air drying unit
- Magnetic strip for automatic load carrier recognition (excl. ML magnets)
- H 500, W 530, D 645 mm

#### Supplied as standard with:

- 9 spring-loaded 12 mm supports for breathing tubes, max. 1.5 m
- 3 x E 433 holders for silicone breathing tubes
- 6 x E 466 injector nozzles for breathing bags
- 15 x E 496 injector nozzles, e.g. for Guedel or endotracheal tubes
- 1 magnetic strip for automatic load carrier recognition (excl. magnets)
- 3 x E 495 injector nozzles
- 1 x E 431 injector nozzle for breathing bag
- 1 x E 430 1/3 mesh tray for sundry small items
- 1 x E 434 holder for paediatric breathing
- 1 x E 432 holder for corrugated breathing tubes

Mat. no. 7 026 590 Art. no. 69.5515.02



#### E 497 adapter

- Nozzle, small
- For double- and single-tubed laryngeal masks
- H 195 mm, with side pipe, 109 mm
- Illustration on right Mat. no. 7 759 470

Art. no. 69.7497.01

#### E 498 adapter

- Nozzle, large
- For double- and single-tubed laryngeal masks
- H 238 mm, with side pipe, 120 mm
- Illustration on left

Mat. no. 7 760 120 Art. no. 69.7498.01



#### E 466 accessory for AN load carrier

Injector nozzle for breathing bag, 8 x 333 mm

Mat. no. 4 966 620 Art. no. 69.7466.01

# E 496 accessory for AN load carrier (no illustration)

Injector nozzle for intubation material,
 4 x 120 mm

Mat. no. 5 235 490 Art. no. 69.7496.01



# Washer-disinfectors G 7825 und G 7826

In terms of size and capacity, G 7825 and G 7826 washer-disinfectors represent the mid-size range among the three Miele series catering for CSSD applications. The characteristic feature of this series is the high throughput capacity of up to 10 DIN mesh trays per cycle and an installation footprint of only **900 mm (width)** x 750 mm (depth).



Illustration shows E 701/2 load carrier for 10 mesh trays with the E 702 top module

| Performance comparison<br>Washer-disinfectors | Unit width/depth | Doors                 | Cabinet Useable dimensions H/W/D Useable volume | Approx.<br>capacity per cycle                |
|---|------------------|-----------------------|---|--|
| PG 8527 / PG 8528                             | 1150/870 mm      | Vertical-rise door(s) | 675/650/800 mm<br>351 I                         | 18 DIN mesh trays<br>7 AN sets<br>3 MI sets  |
| G 7823 / G 7824                               | 900/770 mm       | Bottom-hinged door(s) | 510/530/620 mm<br>168 l                         | 8 DIN mesh trays<br>3 AN sets:<br>2 MI sets  |
| G 7825 / G 7826                               | 900/750 mm       | Bottom-hinged door(s) | 683/541/610 mm<br>225 l                         | 10 DIN mesh trays<br>4 AN sets:<br>2 MI sets |

### Machine features and specifications



#### **Versions**

- G 7825: Front-loading unit with single bottom-hinged door
- G 7826: Two-door barrier washerdisinfector for installation in a diaphragm wall

#### Capacity per cycle

10 DIN mesh trays
 or 4 DIN containers 150 x 300 x 600 mm
 or 2 DIN containers, 300 x 300 x 600 mm
 or 4 AN sets
 or 2 MI sets

#### Design

- Stand-alone or side-by-side installation
- Width 900 mm
- Modular approach with customised features to meet individual requirements
- Single-chamber system for washing, disinfection and drying
- Service-friendly design
- Low heat and sound emissions thanks to double insulation

#### Cleaning technology

- Hygienic freshwater system with fresh water intake in each programme stage
- Cleaning, disinfection and drying in a closed, single-cabinet system
- 2 spray arms in cabinet for thorough cleaning of instrument surfaces
- Spray arms with high water impact force
- Full water jet access, ensuring optimum results
- Thorough cleaning of lumens thanks to injector system
- Direct docking of load carriers to water circuit

#### Standard features

- 2 powerful circulation pumps
- Triple filtration with large surface filter, coarse filter and micro-fine filter
- Filters in inlet hoses
- Flowmeter to monitor water intake quantities
- 1 dump valve

#### **Dispenser systems**

- 2 dispenser pumps for liquid detergent and acidic agents
- Extendible (further 2 internal dispenser pumps retrofittable)



#### Water softeners

• High-capacity water softener (external)

#### Drainage

• 2 drain pumps

#### **Controls**

- Freely programmable controls PROFITRONIC
- 64 programme slots
- 17 standard washing and disinfection programmes
- 8 service programmes
- 39 vacant programme slots
- User interface with local-language display
- Display of programme selection and programming dialogs, countdown time, faults and operating hours.
- Compilation of new programmes using machine controls or using PC/laptop via optical interface

#### Interfaces

- Serial RS 232 interface for process documentation and remote service
- Optical interface for service and maintenance

#### **Printer**

• Integrated 6-pin printer to record key process parameters


#### Safety features

- Electric door lock
- Programme recontinuation in event of power outage
- Peak-load negotiation
- Optical and acoustic signal at end of programme
- 2 separate sensors for temperature monitoring and control
- Port for simple positioning of sensors in the wash cabinet for process validation
- Sensors in cabinet and a magnetic strip on load carriers for automatic assignment of programmes to load carriers

# Modular machine concept Optional extras







#### SBW/2 Plinth/floor tray

- Frame with integrated stainless-steel floor tray
- Cut-outs for steam and plumbing connections, dump valve, dual drain valve and floor tray drain
- A single floor tray can be provided on site for machines installed in a row.
- Plinth fascia (front and rear) with 8 mm recess, flush at sides
- H 100, W 900, D 734 mm

### SBWR/2 Roller plinth/floor tray for G

- Roller plinth/floor tray
- Castors allow machine to be pulled forward for servicing
- For G 7825 with drain pump
- Utility supply access through ceiling
- H 100, W 900, D 734 mm

# PRT/1 Printer for process documentation

- 6-needle printer and RS 232 serial interface
- Paper width: 58 mm
- For installation on infeed side of G 7825 or outfeed side of G 7826
- Recording of following parameters during programme:
   Date and Mach. no., Prog. no., prog. name, starting and finishing time,
- dispensing concentration, dispensing temperature and pumps 1–4, target temperature reached (wash/dry) with times, all faults (e.g. 'water inlet'), all incidents of manual intervention (start, stop, power outage)
- Meets process parameters with respect to temperature and holding time

# Machine versions, components, accessories

#### DK 25/26

#### Steam condenser/heat exchanger

- For G 7825/G 7826 with water cooling (only dehumidified air should be vented to air conditioning system)
- No water consumption when connected to on-site cooling water circuit, max. water pressure 8 bar
- Or cold water supply (water consumption)
- On-site installation
- Reduction of air discharge temperature to approx. 30-35°C
- Reduction in relative humidity to approx.
  60-70%

#### **MAV 25/26**

# Installation kit/Top-box panelling for steam condenser

- Service hatches, with lock, for both sides of diaphragm wall, stainless steel
- Ventilation grille on the infeed side
- Diaphragm wall panelling between the top of the top-box panelling and the ceiling must be provided on site
- H 430, W 900, D 750 mm

#### MVA Installation kit/Floor anchors

- 4 feet, fittings, plugs
- Required when machine is installed without SBW/2

#### OXI-NRS1 - G 7825/G 7826 OxiVario

OxiVario retrofit kit for machines from Serial no. 74305624 incl. on-site installation by Miele Service.

#### OXI-NRS2 - G 7825/G 7826 OxiVario

OxiVario retrofit kit for machines from Mach. no. 53078029 to Mach. no. 74305623, incl. OXI-NRS1 - G 7825/G 7826 and on-site installation by Miele Service.

#### DOS 10/30

#### Integrated dispenser pump

- For surfactant or neutralising agent, 10 ml/30 secs.
- Complete with hoses and siphon (330 mm) for 10 I canisters

#### DOS 60/30

#### Integrated dispenser pump

- For liquid disinfectant or detergent, 60 ml/30 secs.
- Complete with hoses and siphon (330 mm) for 10 I canisters

| G 7825                                | Features                                      | Mat. no. | Art. no.   |
|---------------------------------------|---|----------|------------|
| Electric                              | AE PT EL AV                                   | 5267200  | 62.7825.20 |
|                                       | AE TA PT EL AV                                | 5277520  | 62.7825.21 |
|                                       | AE TA BO PT EL AV                             | 5277530  | 62.7825.22 |
|                                       | AE TA PT EL AP                                | 5544310  | 62.7825.24 |
|                                       | AE TA GS PT EL AP                             | 5769820  | 62.7825.27 |
|                                       | AE TA BO PT EL AP                             | 5430590  | 62.7825.25 |
| Steam                                 | AE BO PT D AV                                 | 5267170  | 62.7825.10 |
| Steam                                 | AE TA BO PT D AP                              | 5430620  | 62.7825.10 |
| G 7826                                | AE IA BO PI D AP                              | 5430620  | 02.7023.13 |
|                                       | AF TA DT FL AV                                | F077000  | 00 7000 01 |
| Electric                              | AE TA CO DT EL AV                             | 5277600  | 62.7826.21 |
|                                       | AE TA BO DT EL AV                             | 5769360  | 62.7826.18 |
|                                       | AE TA BO OF EL AV                             | 5277610  | 62.7826.22 |
|                                       | AE TA BO GS PT EL AV                          | 5769370  | 62.7826.19 |
|                                       | AE TA PT EL AP                                | 5331670  | 62.7826.24 |
|                                       | AE TA BO PT EL AP                             | 5331680  | 62.7826.25 |
|                                       | AE TA BO GS PT EL AP                          | 5769410  | 62.7826.28 |
| Steam/Electric                        | AE BO PT D/EL AV                              | 5267250  | 62.7826.30 |
|                                       | AE TA BO PT D/EL AV                           | 5277630  | 62.7826.32 |
|                                       | AE TA BO GS PT D/EL AV                        | 5769640  | 62.7826.37 |
|                                       | AE TA BO PT D/EL AP                           | 5331700  | 62.7826.35 |
| OxiVario                              | AE TA BO OXI PT EL AV                         | 6407800  | 62.7826.44 |
| CAIVANG                               | TE IT DO ONIT I LETY                          | 0101000  | 02.7020.11 |
| Accessories                           | Features                                      | Mat. no. | Art. no.   |
| SBW/2                                 | Plinth/drip tray for G 7825/G 7826            | 5238130  | 69.3710.02 |
| SBWR/2                                | Plinth/drip tray, on castors for G 7825       | 5653130  | 69.3710.04 |
| MVA                                   | Anchoring kit for G 7825/G 7826               | 5318010  | 69.2100.04 |
| DK 25/26                              | Heat-exchanger steam condenser                | 6600620  | 69.2300.05 |
| MAV 25/26                             | Top-box panelling kit for G 7825/G 7826       | 6600590  | 69.2100.08 |
| OXI-NRS1-G 7825/G 7826                | OxiVario retrofit kit, including installation | 6444680  | 69.2310.03 |
| OXI-NRS2-G 7825/G 7826                | OxiVario retrofit kit, including installation | 6458960  |            |
| PRT/1                                 | Printer for process documentation             | 5400800  | 69.2211.02 |
| DOS 10/30                             | Retrofittable dispenser pump                  | 5267410  | 69.2250.02 |
| DOS 60/30                             | Retrofittable dispenser pump                  | 5267420  | 69.2250.03 |
| Cabinet – ASI 316L (1.4404) grade sta |   |          |            |

# Technical data

| 5 (1 door) G 7826 (2 door)  -  -  -  -  -  -  -  -  -  -  -  -  - | with bottom-hinged door bottom-hinged doors tallations I, max. temperature 93°C docking for cleaning and drying of lumened instruments |
|---|--|
| •<br>•<br>•<br>00° 300/400*                                       | bottom-hinged doors tallations In max. temperature 93°C docking for cleaning and drying of lumened instruments                         |
| •<br>•<br>00* 300/400*<br>•                                       | n, max. temperature 93°C<br>docking for cleaning and drying of lumened instruments   |
| •<br>•<br>•<br>•<br>•   | n, max. temperature 93°C<br>docking for cleaning and drying of lumened instruments   |
| 00* 300/400*  |  |
| •   | FO   |
| •   | s [Qmax. I/min]  |
| •   | nes  |
|   | rogrammable  |
| •   | ds   |
|   |  |
| •   | tion   |
| •   | process documentation  |
| •   | automatic load carrier recognition   |
| •   | abled  |
|   | is   |
| •   | 0 bar flow pressure (200-1000 kPa) (max. 4°dH)   |
| •   | 0 bar flow pressure (200–1000 kPa) (max. 4°dH)   |
| •   | water, 2–10 bar flow pressure (200–1000 kPa)   |
| •   | ith 34" threaded union   |
| •   | odour trap to be fitted on site  |
| 0   | 22, odour trap to be provided on site  |
|   | tion: Electric heating   |
| •   | Iz   |
| 9.0   |  |
| 7* 0.3/0.7*   | kW]  |
| 10.0  | o TA [kW]  |
| 10.0  | th electric TA [kW]  |
| 3 x 16  | an olocatio in the tri   |
| 0 X 10  | tion: Steam heating  |
| •   | Iz   |
| 7* 0.3/0.7*   | kW]  |
| 1.65  | o TA [kW]  |
| 9.0   | th electric TA [kW]  |
| 3 x 16  |  |
| •   | G ½" (DN 15)   |
| •   | e 250–1000 kPa on models with TA/E   |
| •   | e 600-800 kPa on models with TA/D  |
| •   | onnection 600–1200 kPa   |
|   | tion: Steam/Electric (convertible)   |
| •   | Z  |
| 9.0   | ···  |
| 7 0.3/0.7   | kW]  |
| 10.0  | th electric TA [kW]  |
| 3 x 16  |  |
| •   | G ½" (DN 15)   |
| •   |  |
| •   |  |
| •   |  |
| •   | onnection 600-1200 kPa   |
|   | G ½" (DN 15) e 250–1000 kPa on models with TA/E e 600-800 kPa on models with TA/D onnection 600–1200 kPa                               |

# Technical data

| Nasher-disinfectors   | G 7825 (1 door) | G 7826 (2 door) |
|---|-----------------|-----------------|
| Dispenser systems   |                 |                 |
| x DOS 10/30 dispenser pump for liquid acidic agents   | •               | •               |
| x DOS 60/30 dispenser pump for liquid detergent   | •               | •               |
| x 10 I supply canisters   | •               | •               |
| Space for 3 x 5 I supply canisters  | •               | •               |
| Connection options (retrofittable by Service)   |                 |                 |
| OOS 10/30 dispenser pump for neutralising agent   | 0               | 0               |
| OOS 60/30 dispenser pump for detergent/chemical disinfectant                                  | 0               | 0               |
| Orying unit/side-channel compressor   |                 |                 |
| an [kW]   | 2 x 0.9         | 2 x 0.9         |
| leater bank, depending on model [kW]  | 2 x 3.6         | 2 x 3.6         |
| otal rated load, depending on model [kW]  | 9               | 9               |
| Air throughput [m³/h]   | 250             | 250             |
| Temperature selection in 1° increments  | 60–115          | 60–115          |
| Fime selection in 1-minute increments   | 1–240           | 1–240           |
| 2 x pre-filters EU 4, filter rating > 95%, filter life 200 h                                  | •               | •               |
| I particulate/HEPA filters H 13, filter rating > 99.95% (DIN EN 1822-1998), filter life 500 h | •               | •               |
| Dimensions/Weight   |                 |                 |
| External dimensions H incl. plinth tray [mm]  | 1974            | 1974            |
| external dimensions H incl. plinth tray and top-box panelling [mm]                            | 2404            | 2404            |
| External dimensions M/D [mm]  | 900/750         | 900/750         |
| Useable cabinet dimensions H/W/D [mm]   | 683/541/610     | 683/541/610     |
|   |                 | ,,              |
| Overall cabinet dimensions H/W/D [mm]   | 900/567/610     | 900/567/610     |
| Oocking height above floor (with plinth) [mm]   | 850             | 850             |
| Weight [kg]   | 360             | 360             |
| Casing  |                 |                 |
| Stainless steel (AE) Fest certificates  | •               | •               |
|   |                 |                 |
|   |                 |                 |
| e = Standard, o = Optional  |                 |                 |

### Load carrier with 2 - 3 levels



E 775/1 load carrier (empty)

- For inserts on 2 levels
- Built-in spray arm
- Vertical clearances from bottom: Level 1: H 304, W 482, D 590 mm Level 2: H 290, W 488, D 546 mm
- Magnetic strip for automatic load carrier recognition (excl. ML magnets)
- H 400, W 530, D 600 mm

Mat. no. 7 765 730 Art. no. 69.5775.03



#### E 735/2 load carrier (empty)

- For inserts on 3 levels
- 2 built-in spray arms
- Clearances from bottom:
   Level 1: H 203, W 482, D 590 mm
   Level 2: H 203, W 488, D 546 mm
   Level 3: H 133, W 488, D 546 mm
- Magnetic strip for automatic load carrier recognition (excl. ML magnets)
- H 552, W 530, D 600 mm

Mat. no. 7 765 710 Art. no. 69.5735.03 Illustration shows sample load with E 135 containers for baby bottles and E 364 inserts for teats.


E 750 insert frame required for E 135 container for baby bottles. (cf. Page 17)

#### E 775/2 load carrier

- Load capacity:
- 4 DIN mesh trays, 485 x 254 x 50 mm or 2 DIN mesh trays, 540 x 254 x 50 mm or 2 ISO mesh trays, 485 x 344 x 50 mm
- E 750 insert frame required for E 135 container for baby bottles.

#### E 735/2 load carrier

- Load capacity:
- 6 DIN mesh trays, 485 x 254 x 50 mm or 3 DIN mesh trays, 540 x 254 x 50 mm or 3 ISO mesh trays, 485 x 344 x 50 mm



#### A 7/1 insert

- Perforated plate
- 7 x 7 mm perforations
- 3 mm ridge
- For E 775/1, E 735/1 and E 701/1
- H 1, W 543, D 473 mm

Mat. no. 6 097 000 Art. no. 69.5007.02

### Load carrier with 4 - 5 levels





#### E 702 top module

- Level 5 for E 701/2 load carrier
- For 2 DIN mesh trays
- H 160, W 530, D 560 mm

Mat. no. 5 221 490 Art. no. 69.5702.01

# E 701/2 load carrier with drying connection (empty)

- With 4 levels for DIN/ISO mesh trays
- 3 built-in spray arms
- Clearances from bottom:
  Level 1: H 87, W 482, D 590 mm
  Level 2/3: H 87, W 488, D 546 mm
  Level 4 (clearance without Level 5 with E 702): H 223, W 488, D 546 mm
  (Clearance with Level 5, cf. below)

Mat. no. 7 765 700 Art. no. 69.5701.03

# E 701/2 load carrier with E 702 (illustrated)

Level 4: H 87, W 488, D 546 Level 5: H 81, W 488, D 546

- Load capacity:
  8 DIN mesh trays, 485 x 254 x 50 mm
- 4 DIN mesh trays, 540 x 254 x 50 mm or
- 4 ISO mesh trays, 485 x 344 x 50 mm
- Connection for hot-air drying unit
- Magnetic strip for automatic load carrier recognition (excl. ML magnets)
- H 461, W 530, D 600 mm

### Exclusive to Miele

• Max. load per mesh tray 10 kg





- For 7 lumened instruments
- H 170, W 80, D 555 mm

Mat. no. 7 025 350 Art. no. 69.5703.02



#### E 703/1 includes:

- 3 x E 362 blanking screws
- 2 x E 447 female adapter, for male LuerLock
- 2 x E 449 male adapter, for female LuerLock
- 4 x injection nozzles 4 mm, length 110 mm (Mat. no. 4175030)
- 3 x injection nozzles 4 mm, length 120 mm with spring support (Mat. no. 4398951)
- 3 x E 980 angled adapters

### Inserts for shoes, kidney dishes, bowls



#### E 730 insert for 20 theatre shoes

- For E 550 and E 775/2 load carriers
- 20 holders, 280 mm
- For up to 20 pairs of theatre shoes, max. Size 44
- H 290, W 470, D 540 mm

Mat. no. 5 271 020 Art. no. 69.5730.01



#### E 484 insert 1/1 (empty)

- For various utensils
- 1.4 mm wire gauge
- Mesh size: 8 mm
- For holders
- 4 x E 485 for 9 kidney dishes or
- 4 x E 486 for 4 bowls or
- 4 x E 487 for 16 theatre shoes

Oi

- 11 x E 489 universal holders e.g. for insoles
- H 65 (150), W 470, D 480 mm

Mat. no. 5 109 190 Art. no. 69.5484.01



#### E 484 insert with E 487

- E 484 can be fitted with 4 x E 487 holders,
- each for 4 theatre shoes, max. Size 44
- Dimensions E 487 H 280, W 464, D 10 mm

#### E 487

Mat. no. 5 109 250 Art. no. 69.5487.01



#### E 484 insert with E 489

- E 484 can be fitted with 11 x E 489 universal holders, e.g. for insoles
- Dimensions E 489:
   H 60, W 464, D 10 mm

#### E 489

Mat. no. 5 337 120 Art. no. 69.5489.01



#### E 484 insert with E 485

• E 484 can be fitted with 4 holders E 485 for 12 kidney dishes

#### E 485

Mat. no. 5 109 210 Art. no. 69.5485.01



#### E 484 insert with E 486

 E 484 can be fitted with 4 holders E 486 for 4 bowls

#### E 486

Mat. no. 5 109 230 Art. no. 69.5486.01

### Inserts and load carriers for containers



E 484 insert with E 488

• E 484 can be fitted with E 488 for breathing masks

#### E 488

Mat. no. 5 109 260 Art. no. 69.5488.01



E 936 insert for boom lamp handles

- For various lamp handles
- 15 receptacles
- L 500, W 245, H 90 mm

Mat. no. 7 759 480 Art. no. 69.5936.01

#### Supplied as standard with E 713:

- 2 x E 447 female adapters, for male LuerLock
- 3 x injection nozzles, 3 mm, Length 24 mm (Mat. no. 2829560)

Mat. no. 5 221 550 Art. no. 69.5713.01



E 710/1 load carrier (empty)

- For 4 DIN containers and lids
- Container size: 150 x 300 x 600 mm
- Lid size: 30 x 300 x 600 mm
- Magnetic strip for automatic load carrier recognition (excl. ML magnets)
- H 450, W 530, D 600 mm

Mat. no. 5 848 260 Art. no. 69.5710.02



E 711/1 load carrier (empty)

- For 2 DIN containers and lids
- Container size: 300 x 300 x 600 mm
- Lid size: 30 x 300 x 600 mm
- Magnetic strip for automatic load carrier recognition (excl. ML magnets)
- H 455, W 530, D 600 mm


Mat. no. 7 029 790 Art. no. 69.5711.02



E 713 load carrier with drying connector (empty)

- For 2 DIN containers,2 lids and 4 mesh trays on 3 levels
- 2 built-in spray arms
- Clearances from bottom:
   Container size: 150 x 300 x 600 mm
   Lid size: 30 x 300 x 600 mm
   Level 2: H 120, W 488, D 546 mm
   Level 3: H 112, W 488, D 546 mm
- Connections for 5 lumened instruments
- Connection for hot-air drying unit
- Magnetic strip for automatic load carrier recognition (excl. ML magnets)
- H 575, W 530, D 600 mm

### Load carriers for MI instruments





#### E 499 adapter

- Adapter for suction and irrigation cannulae with various connectors, and Verres cannulae
- Compatible with carriages with injector nozzles

Mat. no. 7 759 500 Art. no. 69.5499.01

### E 705/2 mobile injector unit with drying connector

- For MI instruments max. length 600 mm
- Accommodates up to 42 lumened instruments
   (22 connections to top and 20 sloping connections)
- Clearances from bottom:
   Level 1: H 105, W 482, D 590 mm
   Level 2: H 75, W 488, D 546 mm
   Level 3: H 365, W 488, D 546 mm
- Connection for hot-air drying unit
- Magnetic strip for automatic load carrier recognition (excl. ML magnets)
- H 670, W 530, D 600 mm

#### Supplied as standard with:

- 3 x E 336 MIBO irrigation sleeves for MI instruments, 11 mm, Length 121 mm
- 2 x E 362 blanking screws
- 15 x E 442 irrigation sleeves for MI instruments, 4–8 mm
- 5 x E 443 irrigation sleeves for MI instruments, 8–12 mm
- 1 x E 444 insert/spiral rack for fibre optic cables or tubing
- 1 x E 445, 12 caps for irrigation sleeves/MI instruments, 6 mm
- 1 x E 446, 12 caps for irrigation sleeves/MI instruments, 10 mm
- 3 x E 447 female adapters, for male LuerLock
- 6 x E 448 silicone hoses, 300 mm
- 5 x E 449 male adapter, for female LuerLock
- 1 x E 451 insert 1/6 mesh tray, H 55, W 150, D 225 mm

- 3 x E 452 injector nozzles, 2.5 x 60 mm
- 8 x E 453 injector nozzles, with spring support, 4 x 100 mm
- 6 x E 454 injector nozzles for trocar sleeves, 10-15 mm
- 4 x E 456 spring stays for MI instruments
- 3 x E 464 adapter for E 454 injector nozzle for trocar sleeves, 20 mm
- 2 x E 472 clamp for injector nozzle, 4.0 mm
- 2 m silicone hose, 5 mm Part no. 4822890

Mat. no. 7 029 810 Art. no. 69.5705.03

# Load carrier for MI instruments and reprocessing dental transmission instruments



## E 774/1 mobile injector unit with drying connector

- Modular system for 2-3 MI OP sets or dental transmission instruments
- For lumened instruments in 4 modules with integrated nozzles/adapters (not included)
- Applications:
- E 903/1 module for MI instruments/ urology
- E 905/1 module for short MI instruments E 906/1 module for long MI instruments 3  $\times$  E 444 drums for lightguides and irrigation hoses
- 1 x E 460 insert for rigid fibre optics 1 x E 457 insert for separable MI instruments or
- 1 x E 142 DIN mesh tray
- Connection for hot-air drying unit
- Magnetic strip for automatic load carrier recognition (excl. ML magnets)
- H 605, W 530, D 600 mm

#### Supplied as standard with:

- 4 x E 447 female adapters, for male LuerLock
- 6 x E 362 blanking screws

Mat. no. 7 040 190 Art. no. 69.5774.02

### Load carriers for micro-instruments



Geuder AG recommends that their current range of instruments be reprocessed in a Miele washer-disinfector.





#### E 729 load carrier with drying connector

- For 4–5 instrument sets (cataracts)
- Load carrier with 3 levels for microinstruments and accessories
- 2 levels for 4 DIN mesh trays
- 1 level with 16 injector nozzles (with hose and male LuerLock) and 30 injector nozzles (15 male and 15 female LuerLock adapters)
- Connection for E 478 on feed pipe
- Vertical clearance (bottom to top):

Level 1: 124 mm Level 2: 114.5 mm Level 3: 160 mm

- Connection for hot-air drying unit
- Magnetic strip for automatic load carrier recognition (excl. ML magnets)
- H 471, W 530, D 600 mm

#### Supplied as standard with:

- 1 x E 476 adapter
- 1 x E 477 stopper
- 8 connectors LuerLock female/female (E 790)
- 8 connectors LuerLock female/male (E 791)

Mat. no. 6 517 580 Art. no. 69.5729.01

#### E 478/1 adapter

• For 4 narrow-lumened cannulae (Sautter cannulae)

Mat. no. 7 292 920 Art. no. 69.7478.02



#### FP ED

- Stainless-steel filter insert for E 478/1
- Ø 30 mm

Mat. no. 7 290 770 Art. no. 69.7512.02

### Load carriers for anaesthetic equipment



### E 715/1 mobile injector unit with drying connector

- For approx. 4 AN sets or
   12 breathing tubes, max. length 1.5 m,
   and various accessories
- Tubes wound around spiral drum
- Connection for hot-air drying unit
- Magnetic strip for automatic load carrier recognition (excl. ML magnets)
- H 645, W 530, D 600 mm

#### E 715/1 supplied as standard with:

- 1 x E 430 insert 1/3 mesh tray, H 40, W 150, D 445 mm (supplied)
- 1 x E 431 injector nozzle for breathing bags,
- 8 x 193 mm (supplied)
- 1 x E 432 holder for 3 or 4 corrugated breathing tubes (supplied)
- 3 x E 433 holders for 3 or 4 silicone breathing tubes
- 1 x E 434 holder for 3 or 4 paediatric breathing tubes (supplied)
- 6 x E 466 injector nozzles for breathing bags,
- 8 x 333 mm
- 15 x E 496 injector nozzles, 4 x 120 mm

Mat. no. 6 491 890 Art. no. 69.5715.02



#### E 497 adapter

- Nozzle, small
- For double- and single-tubed laryngeal masks
- H 195 mm, with side pipe, 109 mm
- Illustration on right Mat. no. 7 759 470 Art. no. 69.7497.01

#### E 498 adapter

- Nozzle, large
- For double- and single-tubed laryngeal masks
- H 238 mm, with side pipe, 120 mm
- Illustration on left Mat. no. 7 760 120

Art. no. 69.7498.01



#### E 466 accessory for AN load carrier

Injector nozzle for breathing bag, 8 x 333 mm

Mat. no. 4 966 620 Art. no. 69.7466.01

# E 496 accessory for AN load carrier (no illustration)

Injector nozzle for intubation material,
 4 x 120 mm

Mat. no. 5 235 490 Art. no. 69.7496.01

# Article overview

| Article             | Page       |
|---------------------|------------|
| A 6                 | 39         |
| A 7/1               | 68         |
| A 16                | 42         |
| A 18                | 42         |
| A 502               | 36         |
| A 600               | 36         |
| A 602               | 36         |
| A 805               | 36         |
| A 806               | 36         |
| A 807               | 36         |
| Cover plate (23/24) | 51         |
| DE 23               | 51         |
| DK 23/24            | 50         |
| DK 25/26            | 65         |
| DK 27/28            | 21         |
| DOS 10/30 (23/24)   | 51         |
| DOS 10/30 (25/26)   | 65         |
| E 131/1             | 41         |
| E 142               | 39, 42     |
| E 143               | 42         |
| E 146               | 42         |
| E 336               | 38         |
| E 362               |            |
| E 363               | 43         |
|                     |            |
| E 373               | 40, 41     |
| E 374               | 40         |
| E 378               | 42         |
| E 379               | 42         |
| E 416               | 40         |
| E 417/1             | 40         |
| E 430/1             | 42         |
| E 441/1             | 41         |
| E 442               | 38         |
| E 443               | 38         |
| E 444               | 43         |
| E 445               | 38         |
| E 446               | 38         |
| E 447               | 38         |
| E 448               | 38         |
| E 449               | 38         |
| E 451               | 41         |
| E 452               | 38         |
| E 453               | 38         |
| E 454               | 38         |
| E 456               | 38         |
| E 457               | 43         |
| E 460               | 43         |
| E 464               | 38         |
| E 466               | 31, 59, 75 |
| E 467               | 38         |
| E 468               | 42         |
| E 469               | 38         |
| E 471               | 38         |
| E 472               | 38         |
| E 473/1             | 41, 43     |
| E 475               | 38         |
| E 476               | 39         |
| E 477               | 39         |
| E 478/1             | 30, 58, 74 |

| Article | Page      | е |
|---------|-----------|---|
| E 479   | 39        |   |
| E 484   | 54, 70, 7 | 1 |
| E 485   | 70        | 0 |
| E 486   | 70        | 0 |
| E 487   | 54, 70    | 0 |
| E 488   | 7         | 1 |
| E 489   | 70        | 0 |
| E 492   | 2         | 7 |
| E 496   | 31, 59, 7 | 5 |
| E 497   | 31, 59, 7 | 5 |
| E 498   | 31, 59, 7 | 5 |
| E 499   | 39, 50    | 6 |
| E 515/1 | 1 59      | 9 |
| E 517   | 38        | 8 |
| E 524/2 | 2 5       | 5 |
| E 525/2 | 2 5       | 1 |
| E 526/1 | 1 TA 50   | 6 |
| E 527/2 | 2 5       | 5 |
| E 528/2 | 2 5       | 5 |
| E 529   | 58        | 8 |
| E 550   | 54        | 4 |
| E 555/2 | 2 5       | 5 |
| E 574   | 5         | 7 |
| E 701/2 | 2 TA 69   | 9 |
| E 702   | 54, 69    | 9 |
| E 703/1 | 1 55, 69  | 9 |
| E 705/2 | 2 TA 72   | 2 |
| E 710/1 | 1 7       | 1 |
| E 711/1 | 1 7       | 1 |
| E 713 T | ΓA 7      | 1 |
| E 715/1 | 1 TA 75   | 5 |
| E 729 T | ΓA 7-     | 4 |
| E 730   | 54, 60    | 0 |
| E 735/2 | 2 6       | 8 |
| E 750   | 5         | 5 |
| E 774/1 | 1 TA 73   | 3 |
| E 775/1 | 1 6       | 8 |
| E 790   | 39        | 9 |
| E 791   | 39        |   |
| E 792   | 39        | 9 |
| E 803   | 40        | 0 |
| E 900-4 | 4/2 24    | 4 |
| E 900-5 |           |   |
| E 901/2 |           | 5 |
| E 902/1 |           |   |
| E 902/3 |           |   |
| E 903/1 |           |   |
| E 905/1 | 1 3-      | 4 |
| E 906/1 |           |   |
| E 907/1 |           | 3 |
| E 908/1 |           |   |
| E 909 T |           |   |
| E 910/4 |           |   |
| E 911/3 |           |   |
| E 912/2 |           |   |
| E 913/3 |           |   |
| E 915/3 |           |   |
| E 919   | 33        |   |
| E 926   | 33        |   |
| E 929/1 | 1 TA 30   | 0 |
|         |           |   |

| Article                   | Page       |
|---------------------------|------------|
| E 930                     | 26         |
| E 931                     | 26         |
| E 934/1                   | 28         |
| E 935/2                   | 24         |
| E 936                     | 27, 55, 71 |
| E 941 TA                  | 29         |
| E 948                     | 29         |
| E 949                     | 29         |
| E 955                     | 29         |
| E 975/2                   | 24         |
| E 980                     | 25         |
| E 984                     | 26         |
| E 985                     | 27         |
| E 986                     | 27         |
| E 987                     | 26         |
| E 989                     | 26         |
| FP ED                     | 30, 58, 74 |
| MAV 23/24                 | 51         |
| MAV 25/26                 | 65         |
| MAV 27/28                 | 21         |
| MF/3                      | 44         |
| MF 27/28-1                | 44         |
| ML                        | 51         |
| OrthoVario module (23/24) | 51         |
| OxiVario module (23/24)   | 51         |
| MVA                       | 27         |
| ORTHO (23/24)             | 51         |
| ORTHOVARIO (27/28)        | 21         |
| OXIVARIO                  | 21         |
| OXI-NRS1-G 7823/G 7824    | 51         |
| OXI-NRS1-G 7825/G 7826    | 65         |
| OXI-NRS2-G 7825/G 7826    | 65         |
| Printer                   | 20         |
| PRT/1                     | 51, 64     |
| SBW                       | 20         |
| SBW/1                     | 20         |
| SBW/2                     | 64         |
| SBW/4                     | 50         |
| SBWR                      | 20         |
| SBWR/2                    | 64         |
| SBWR/3                    | 50         |
| TA/D                      | 21         |
| TA/E                      | 21         |
| TA-S1                     | 50         |
| TK/1                      | 27         |
| UTS                       | 42         |



# Notes

# Notes



Miele Co. Ltd. Fairacres Marcham Road, Abingdon Oxon OX14 1TW

Tel 0845 330 3618

miele-professional@miele.co.uk

### Miele worldwide:

AUSTRALIA Miele Australia Pty. Ltd. 1 Gilbert Park Drive Knoxfield, VIC, 3180 Tel. +613 97 64 71 30 Telefax +613 97 64 71 49

AUSTRIA Miele Gesellschaft mbH. Mielestraße 1 5071 Wals/Salzburg Tel. +43 50 800-81 481 Telefax +43 50 800-81 429

BELGIUM N.V. Miele S.A. z.5 Mollem 480 Hof te Bollebeeklaan 9 1730 Mollem Tel. +32 2451 14 11 Telefax +32 2451 14 14

CANADA Miele Limited 161 Four Valley Drive Vaughan, ON L4K 4V8 Tel. +1 905 532 2270 Telefax +1 904 532-2290

CHILE
Miele Electrodomesticos Ltda.
Nueva Costanera 4055
Vitacura
Santiago de Chile
Tel. +56 2 95700-00
Telefax +56 2 95700-79

CHINA
Miele (Shanghai) Trading Ltd
Beijing Branch
3F South Tower, Kerry Center
1 Guang Hua Road
Chaoyang District
100020 Beijing
Tel. +96 10 5783 2688
Telefax +86 10 5783 2600

Miele (Shanghai) Trading Ltd. 1-3 Floor, No. 82 Shi Men Yi Road, Jing An District 200040 Shanghai PRC Tel. +86 21 6135 3559 Telefax +86 21 6135 3993

CROATIA
Miele trgovina i servis d.o.o.
Buzinski prilaz 32
10000 Zagreb
Tel. +385/1/66 89 000
Telefax +385/1/66 89 090

CZECH REPUBLIC Miele spol. s.r.o. Holandská 4 63900 Brno Tel. +420 5 4 35 53-1 11 Telefax +420 5 4 35 53-1 19

DENMARK
Miele A/S
2600 Glostrup, Erhvervsvej 2
Tel. +45 43 27 11 00
Tel. Saldy/Professional
+45 43 27 15 00
Tel. Jylland +45 97 12 70 66
Telefax hovednr. +45 43 27 11 09
Telefax Salg/Professional
+45 43 27 15 09

ESTONIA Miele SIA Eesti filiaal A.H. Tammsaare tee 47 11316 Tallinn Tel. +371 2922 3571

FINLAND Miele Oy Porttikaari 6 01200 Vantaa Tel. +3589 87 59 70 Telefax +3589 87 59 72 99

FRANCE Miele S.A.S. 9 avenue Albert Einstein 93151-Le Blanc-Mesnil (Paris) Tel. +33(0)1/49 39 34 44 Telefax +33(0)1/49 39 44 38

GERMANY Miele & Cie. KG Vertriebsgesellschaft Deutschland Carl-Miele-Str. 29 33332 Gütersloh Tel. 0800 22 44 644\* Telefax 0800 33 55 533\* \*Only in Germany

GREAT BRITAIN Miele Co. Ltd. Fairacres Marcham Road, Abingdon Oxon OX14 1TW Tel. 0845 330 3618 Telefax 01235 405636

GREECE Miele Hellas E.P.E. Mesogion 257 154 51 N. Psychiko-Athens Tel. national 01/10 6 79 44 44 international 00 30/2/10 6 79 44 44 Telefax 01/10 6 79 42 09 international 00 30/2/10 6 79 42 09

HUNGARY Miele Kft. Alsó Törökvész út 2. 1022 Budapest Tel. +36 1 8806-400 Telefax +36 1 8006-402

HONG KONG Miele (Hong Kong) Ltd. 41/F - 4101, Manhattan Place 23 Wang Tai Road, Kowloon Bay Tel. +852/2610 1331 Telefax +852/3579 1404

INDIA Miele India Pvt. Ltd. 1st Floor, Copia Corporate Suites, Plot No. 9, Jasola, New Delhi – 110025 Tel. +91 (0)11 469000-00 Telefax +91 (0)11 469000-01

IRELAND Miele Ireland Ltd. 2024 Bianconi Avenue Citywest Business Campus Dublin 24 Tel. +3531 461 07 10 Telefax +3531 461 07-97 ITALY Miele Italia S.r.I. Strada di Circonvallazione, 27 39057 S. Michele-Appiano (BZ) Tel. +39 04 71/66 61 11 Telefax +39 04 71/66 63 50

JAPAN Miele Japan Corp. Meguro Yamate Place 9F 2-10-11 Meguro Meguro-ku, Tokyo 153-0063 Japan Tel. +81(3)5740-0034 Telefax +81(3)5740-0035

KAZAKHSTAN Miele LLP Building 1B (office 602) Al-Farabi Pr. 13 050059 Almaty Tel. +7 727 3111141 Telefax +7 727 3111042

KOREA Miele Korea Limited Miele House 8 fl., 607-10 Yeoksam-dong Gangnam-gu Seoul 135-080 Tel. +82 2 3451 9333 Telefax +82 2 3451 9399

LATVIA Miele SIA Zaubes Str. 9A-38 1013 Riga Tel. +371 66016600 Telefax +371 67889879

LITHUANIA SIA "Miele" Lietuvos filialas Laisvės pr. 125 06118 Vilnius Tel. +371 29223521

LUXEMBOURG Miele s.à r.l. 20, Rue Christophe Plantin 2339 Luxembourg Tel. +352/4 97 11-25 Telefax +352/4 97 11-39

MALAYSIA Miele Sdn Bhd Suite 12-2, Level 12 Menara Kencana Petroleum Solaris Dutamas No. 1, Jalan Dutamas 1 Kuala Lumpur 50480 Tel. +603 6209 0288 Telefax +603 6205 3768

MEXICO Miele S.A. de C.V. Av. Santa Fé 170 German Centre 0-4-2 Col. Lomas de Santa Fé C.P. 01210 México, D.F. Tel. +52/55 85 03 98 70 Telefax +52/55 85 03 98 74

NETHERLANDS Miele Nederland B.V. De Limiet 2 4131 NR Vianen Tel. +31 3 47/37 88 83 Telefax +31 3 47/37 84 29 NEW ZEALAND Miele New Zealand Ltd. Unit L, 10–20 Sylvia Park Road Mt. Wellington, 1060 Auckland Tel. +64 9573 1269 Telefax +64 9573 1268

NORWAY Miele AS Nesbruveien 71 1394 NESBRU Tel. +47/67 17 31 00 Telefax +47/67 17 31 10

POLAND Miele Sp. z. o. o. ul. Gotarda 9 02-683 Warszawa Tel. +48 22/5 48 40 00 Telefax +48 22/5 48 40 10

PORTUGAL Miele Portuguesa, Lda. Av. do Forte, 5 2790-073 Carnaxide Tel. +351/21/42 48-100 Telefax +351/21/42 48-109

ROMANIA Miele Appliances SRL Intrarea Narciselor 8 Otopeni Jud. Ilfov Tel. +40 21 3520777 Telefax +40 21 3520776

RUSSIA 000 Miele CIS Leningradsky Prospekt, 31A, bld. 1 125284 Moscow Tel. +7 495 745-8992 Telefax +7 495 745-8680

SERBIA Miele d.o.o. Beograd Balkanska 2 Beograd – Stari Grad 11000 Tel. +381 11 2055060 Telefax +381 11 2055065

SINGAPORE Miele Pte. Ltd. 163 Penang Road #04-02/03 Winsland House II Singapore 238463 Tel. +65/67 35 11 91 Telefax +65/67 35 11 61

SLOVAKIA Miele s. r. o. Plynárenská 1 82109 Bratislava Tel. +421/2/58 10 31 11 Telefax +421/2/58 10 31 19

SLOVENIA Miele d.o.o. Brnőičeva ulica 41 g 1231 Ljubljana – Črnuče Tel. +386 1563 44-80 Telefax +386 1563 44-90

SOUTH AFRICA Miele (Pty) Ltd. 63 Peter Place Bryanston 2194 Tel. +27(0)11/875-9000 Telefax +27(0)11/875-9035 SPAIN Miele S.A.U. Avda. de Bruselas, 31 28108 Alcobendas (Madrid) Tel. +34/91/6 23 20 00 Telefax +34/91/6 62 02 66

SWEDEN Miele AB Industrivägen 20 Box 1397 171 27 Solna Tel. +46 85 62 29-000 Telefax +46 85 62 29-209

SWITZERLAND Miele AG Limmatstr. 4 8957 Spreitenbach Tel. +4156 4 17 20 00 Telefax +4156 4 17-24 69

TURKEY Miele Elektrikli Aletler Dış. Tic. ve Paz. Ltd. Şti. Güvercin Sokak No. 10 34330 Levent-Istanbul Tel. +90/2 12/3 86-08 00 Telefax +90/2 12/3 25 84 49

UKRAINE Miele LLC Zhylyanskaya Str. 48, 50 A 01033 Kiev Tel. +38 044 4960300 Telefax +38 044 4942285

UNITED ARAB EMIRATES Miele Appliances Ltd. Gold & Diamond Park Sheikh Zayed Road Building No. 6, Offices No. 219 Dubai Tel. +971 4 3418444 Telefax +971 4 3418852

USA Miele, Inc. 9 Independence Way Princeton, NJ 08540 Tel. 001-609/419 43 74 Telefax 001-609/419 18 53

OTHER COUNTRIES: Sales International Miele & Cie. KG Carl-Miele-Straße 29 33332 Gütersloh Tel. +49(0)52 41-89 15 07 Telefax +49(0)52 41-89 5 00